



Cona (VE) Italy, 17/07/2025

TO WHOM IT MAY CONCERN:

**DECLARATION
ABSENCE OF PHTHALATES**

Based on the information currently in our possession, we report what is declared by the Manufacturer of the raw materials used to produce the products listed below:

CL4ARBIPSTS/B Cryogen® Tube 4 MI Self Standing Internal Stella Cap, Graduated, Bar-Code-Sterile-Clearline® In Bulk

Cap (Polyethylene)– Regarding classification of the above product according to the EU Directive 1999/45/EC, which refers to Directive 67/548/EEC and subsequent amendments, reference is made in the SDS/PSIS for the above product. We also certify, that during the manufacturing of the above product, we do not use or intentionally incorporate into it any of the following materials: Plasticisers (e.g. Adipates, ESBO, Phthalates)

Tube (Polypropylene)- The composition of this product meets the relevant requirements of the following documents: Directive 2011/65 and amendments, on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS) which establishes a maximum concentration value up to: **Bis(2-ethylhexyl) phthalate (DEHP) (0,1 %)**

O-Ring (TPV) - With regard to the compliance status of the product referenced above with the regulation(s) identified below the following can be declared:

This product is in compliance with the relevant heavy metals, flame retardants and phthalates requirements of the following regulation:

The concentrations of the following heavy metals (lead, cadmium, mercury & hexavalent chromium) flame retardants [polybrominated biphenyls (PBBs), polybrominated diphenyl ethers (PBDEs)] the following **phthalates [Bis(2-Ethylhexyl) phthalate (DEHP), Benzyl butyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP)]** in this product do not exceed 0.1% by weight for lead, mercury, hexavalent chromium, PBBs, PBDEs & phthalates and 0.01% by weight for cadmium. Traces levels of these substances may be present resulting from the specific characteristics of the raw materials and/or of the manufacturing process.



Biosigma S.p.A.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 | Fax ++39 0426 302228 |
Email info@biosigmaeu.com | www.biosigma.it www.biosigma.com | Cap. Soc. 99.000 € i.v. | C.F. 02173800281 |
Codice SDI: MZO2A0U | P.IVA IT 03328440270 | N. 298577 R.E.A. di Venezia | Registro A.E.E. IT09030000005843 | R.N.P.P.A. IT11080P00002626