











SEZIONE A-A



Max dimensions of
bar-code printing
area 20 x 10 mm.



DETTAGLIO B
SCALA 4 : 1

DIN 16901			SCOSTAMENTI AMMISSIBILI																				
Gruppo Tolleranza	da	0	1	3	6	10	15	22	30	40	53	70	90	120	160	200	250	315	400	500	630	800	
	a	1	3	6	10	15	22	30	40	53	70	90	120	160	200	250	315	400	500	630	800	1000	
	** A) per dimensioni non legate allo stampo										TOLLERANZE GENERALI												** B) per dimensioni legate allo stampo
130	A	± 0.15	± 0.19	± 0.20	± 0.21	± 0.23	± 0.25	± 0.27	± 0.30	± 0.34	± 0.38	± 0.44	± 0.51	± 0.60	± 0.70	± 0.90	± 1.10	± 1.30	± 1.60	± 2.00	± 2.50	± 3.00	
	B	± 0.05	± 0.09	± 0.10	± 0.11	± 0.13	± 0.15	± 0.17	± 0.20	± 0.24	± 0.28	± 0.34	± 0.41	± 0.50	± 0.60	± 0.80	± 1.00	± 1.20	± 1.50	± 1.90	± 2.40	± 2.90	
A	Nuovo fondello a 6 intagli													09-2018		MTE							
Revisione \ Revision	Descrizione \ Description													Data \ Date		Firma \ Signature		Approvato \ Approved					
 Biosigma srl a dominique Dutscher company Via Valletta 6 30010 Cantarana di Cona VE www.biosigma.com tel. 0426 3022241 progettazione@biosigmaeu.com		Materiale \ Material Test-Tube in PP Med Cap in Silicon		Peso \ Weight 8.21		Finitura \ Roughness					Scala \ Scale 2:1		Revisione \ Revision -										
		Codice \ Code CL5ARBIPSTS									Data \ Date		Disegnato da \ Draw by DA										
		Descrizione \ Description		Cryogen tube 5 ml int. thread, self-stand, with printing and bar-code area							Stampo \ Mold		Approvato da \ Approved by DA										