




# flip&grip®

## pipette holder

### Flexible bench or under-shelf storage

Unique pipette holder that gives you the option of using on the bench top, or hanging under shelf holder to help maximize your working space.

- Will hold up to 12 single channel pipette, or a combo of 3 single and 3 multichannel pipettes
- Holds most common brands
- Rubber grips prevent scratching pipette
- Holder rotates for easy pipette selection and storage
- Weighted base
- Fixing screw provided for under shelf option

Part No.	H x W Bench Top Style	H x W Under Shelf Style	Color
<b>HS120026</b>	332 x 193 mm	187 x 193 mm	Blue 
<b>HS120056</b>	332 x 193 mm	187 x 193 mm	Lime 
<b>HS120057</b>	332 x 193 mm	187 x 193 mm	Magenta 



Pipettes can be held in the rubber grips and/or the plunger hooks



Mount base under shelf or counter



Rack rotates for easy pipette selection



Each unit is color coded and labeled with volume range (Both top and side for quick visual confirmation)



### pearl® pipettes adjustable volume

Precision liquid handling with a complete range of single channel pipette volumes to cover all your possible needs.

The sleek ergonomic design and reduced plunger pressure help reduce the chance of Repetitive Stress Injury (RSI), as the pipettes require 2/3rd the operation pressure of most standard pipettes.


Each pipette is tested in accordance with standards and a calibration certificate is supplied.

Manufactured in accordance with ISO08655 and GLP regulations.

- Internal metal parts for increased accuracy and calibration stability
- Easy volume adjustment with locking tab to prevent accidental change
- Separate tip ejector and a large digital volume display
- CE certified
- **Autoclavable**

#### Features:

1. Volume adjustment and low-force dispensing all in one button for convenience
2. Slight downward slant to tip ejector button reduces hand fatigue, incidence of repetitive strain injury (RSI)
3. Isolated volume locking ring prevents accidental volume changes
4. Finger hook
5. 3-position volume display
6. Universal tip cone



## 10 Volume Ranges Available

Part No.	No. Channels	Volume Range	Color
HS120097	Single	0.1–2 µL	Green
HS120098	Single	0.5–10 µL	Lt. Blue
HS120099	Single	2–20 µL	Lt. Purple
HS120100	Single	5–50 µL	Orange
HS120101	Single	10–100 µL	Blue
HS120102	Single	20–200 µL	Yellow
HS120103	Single	30–300 µL	Purple
HS120104	Single	100–1000 µL	Blue
HS120105	Single	500–5000 µL	Pink
HS120106	Single	1000–10000 µL	Red
HS120215	Eight	5–50 µL	Orange
HS120216	Eight	30–300 µL	Purple

Precision metal components for both shaft and barrel ensure reliable performance and long life



8 channel pipette available in 5–50 µL and 30–300 µL



Heathrow Scientific, LLC | 620 Lakeview Parkway, Vernon Hills IL 60061 USA  
Tel.: +1-800-741-4597 | Fax: +1-847-816-5072 | E-mail: info@heatsci.com