

# adapt-a-rack™

## Flexible multi-tube rack

Each slot adapts to accommodate 12 mm–30 mm (5 mL–50 mL) tubes sizes and holds them in a firm upright position

- Made from POM
- Racks can be linked together for additional holding capacity
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Linked racks should be supported underneath when carrying
- Available in solid or a combination of colors for easy coding identification

- **Autoclavable**

Part No.	Tubes	Wells	Wells Dia.	L x W x H	Pkg.	Rack #1 Color	Rack #2 Color
HS120184	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Blue/Green	Green/Blue
HS120185	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Blue/Yellow	Yellow/Blue
HS120186	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Pink/White	White/Pink
HS120187	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Blue/Blue	Blue/Blue
HS120188	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	White/White	White/White
HS120189	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Pink/Pink	Pink/Pink
HS120190	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Yellow/Yellow	Yellow/Yellow
HS120191	5 – 50 mL	4	12 – 30 mm	7.12 x 2.20 x 3.0 in (18.1 x 5.6 x 7.6 cm)	2/pk	Green/Green	Green/Green



Available in solid or a combination of colors




Holds **different size tubes** in each slot



## Universal Tube Rack



Ideal for holding a wide range of tubes

- Made from polypropylene
- Holds up to 12 similar size tubes
- Imprinted references for easy identification
- Side handles for easy transport
- **Not recommended for autoclaving**



Adjustable slot adapts to different diameter tubes

Part No.	Tubes	Wells	Well Dia.	L x W x H	Pkg.	Color
<b>HS24330A</b>	15-50 mL	12	16-30 mm	7.6 x 5.5 x 2 in (19.4 x 14 x 5.1 cm)	1/ea	Blue/Purple <span style="color:blue">■</span> <span style="color:purple">■</span>



Hold and pour feature

## grab rack™

Rack features silicon grips to hold tubes in place

- ABS rack with silicon gripping fingers
- Turn the rack over and the tubes will not fall out
- No assembly required
- **Not recommended for autoclaving**


Part No.	Tubes	Wells	Well Dia.	Dimensions L x W x H	Pkg.	Color
<b>HS24371A</b>	5-10 mL	50	16 mm	10 x 5 x 2 in (25.9 x 12.7 x 4.9 cm)	1/ea	White <input type="checkbox"/>

## Conical Centrifuge

Available in natural and pure black for light sensitive samples

- Attached flat caps for easy one handed use
- 1.5 mL versions have a frosted writing area and a syringe pierceable cap
- 100% virgin polypropylene
- **Disposable**

Certified DNase, RNase and Pyrogen free



Skirted

Large capacity tubes – 15 and 50 mL

- Lot-to-lot traceability with easy tear packaging
- Writing area and flat polyethylene lids for easy labelling
- Imprinted graduations and a self-standing version (only in clear bulk pack)
- Reusable card used for the racked versions, making them more environmentally friendly

	Part No.	Vol.	Dimensions	Graduations	Max RCF Rating	Pkg.	Color
<b>a</b>	<b>HS4426B</b>	15 mL	OD 17 mm Length 120 mm	0.5 mL	11,000 x g	Bulk wrap sterile	Natural <input type="checkbox"/>
	<b>HS4426R</b>	15 mL	OD 17 mm Length 120 mm	0.5 mL	11,000 x g	Racked sterile	Natural <input type="checkbox"/>
<b>b</b>	<b>HS4428</b>	15 mL	OD 17 mm Length 120 mm	0.5 mL	9,400 x g	Bulk wrap sterile	Black <input checked="" type="checkbox"/>
	<b>HS4427B</b>	50 mL	OD 29 mm Length 114 mm	5.0 mL	9,500 x g	Bulk wrap sterile	Natural <input type="checkbox"/>
<b>c</b>	<b>HS4427R</b>	50 mL	OD 29 mm Length 114 mm	5.0 mL	9,500 x g	Racked sterile	Natural <input type="checkbox"/>
	<b>HS4429</b>	50 mL	OD 29 mm + 0.3 mm Length 114 mm	5.0 mL	9,400 x g	Bulk wrap sterile	Black <input checked="" type="checkbox"/>
<b>e</b>	<b>HS4427S</b>	50 mL	OD 29 mm + 0.3 mm Length 114 mm	5.0 mL	9,500 x g	Bulk wrap sterile skirted	Natural <input type="checkbox"/>