

Cardboard Slide Folders

The cardboard slide folders by Heathrow Scientific offer a convenient solution for safely working, transporting and short-term storing of slides.

- Easily organize slides with the index on the front cover that corresponds with the number grid inside
- Recessed slide compartments and inside covers protect specimens. Slides cannot be touched when the cover is closed, reducing contamination
- Front protection flaps are attached with durable tape for repeated use. Fold flaps flat under the slide holder to save space when in use
- Thumb cut in slide compartments enables easy removal of slides with slight pressure at one end
- Compartment holder is made from a double layer of densely pressed cardboard. Cover is a single layer cardboard
- Available in two sizes that hold 20 or 30 standard (26 x 76mm and 1 x 3in) slides
- Disposable


Item No.	Color	Capacity	L x W x H in	L x W x H cm	UOM
HS9902	Black				
HS9913	Blue				
HS9914	Green	20 slides	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9	72/cs
HS120196	Red				
HS120197	Yellow				
HS9930	Black	30 slides	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1	40/cs



Slide Folder Rack

Effortlessly locate, retrieve and keep 20-place slide folders organized

- Made from ABS
- Ideal for histology, microbiology and pathology as well as multiple user research environments
- Rack is lightweight and can be safely transported to multiple work areas
- Rack is portable and refrigerator-safe for cold sample storage
- Multiple folder racks can be stacked vertically to save space in the lab
- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached
- Additional clasps can be purchased



Stay Organized
Includes 5 assorted colored clasps



Stackable

Slide Folder Rack

Item No.	Color	Capacity	L x W x H in	L x W x H cm	UOM
HS99040A	Blue	5 Slide Trays	8.9 x 13.4 x 4.9	22.6 x 34 x 12.5	1/ea

**5 Clasps Included (Clear, Blue, Teal, Purple, Pink)*

Slide Folder Clasps

Item No.	Color	L x W x H in	L x W x H cm	UOM
120653	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	2.2 x 1.7 x 0.6	5.6 x 4.4 x 1.6	20 ea color/bag
120654	Clear <input type="checkbox"/>			100 ea/bag

**Assorted colors (Clear, Gray, Blue, Green, Red)*



120653