

# Cardboard Slide Folders

The cardboard slide folders by Heathrow Scientific offer a convenient solution for safely working, transporting and short-term storing of slides.

- Easily organize slides with the index on the front cover that corresponds with the number grid inside
- Recessed slide compartments and inside covers protect specimens. Slides cannot be touched when the cover is closed, reducing contamination
- Front protection flaps are attached with durable tape for repeated use. Fold flaps flat under the slide holder to save space when in use
- Thumb cut in slide compartments enables easy removal of slides with slight pressure at one end
- Compartment holder is made from a double layer of densely pressed cardboard. Cover is a single layer cardboard
- Available in two sizes that hold 20 or 30 standard (26 x 76mm and 1 x 3in) slides
- Disposable


| Item No. | Color    | Capacity  | L x W x H<br>in  | L x W x H<br>cm   | UOM   |
|----------|----------|-----------|------------------|-------------------|-------|
| HS9902   | Black ■  |           |                  |                   |       |
| HS9913   | Blue ■   |           |                  |                   |       |
| HS9914   | Green ■  | 20 slides | 8.2 x 13.4 x 0.4 | 20.8 x 34.1 x 0.9 | 72/cs |
| HS120196 | Red ■    |           |                  |                   |       |
| HS120197 | Yellow ■ |           |                  |                   |       |
| HS9930   | Black ■  | 30 slides | 8.2 x 12.9 x 0.4 | 28.1 x 32.8 x 1   | 40/cs |



# Slide Folder Rack

Effortlessly locate, retrieve and keep 20-place slide folders organized

- Made from ABS
- Ideal for histology, microbiology and pathology as well as multiple user research environments
- Rack is lightweight and can be safely transported to multiple work areas
- Rack is portable and refrigerator-safe for cold sample storage
- Multiple folder racks can be stacked vertically to save space in the lab
- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached
- Additional clasps can be purchased



**Stay Organized**  
Includes 5 assorted colored clasps



HS99040A



**Stackable**

## Slide Folder Rack

| Item No. | Color | Capacity      | L x W x H<br>in  | L x W x H<br>cm  | UOM  |
|----------|-------|---------------|------------------|------------------|------|
| HS99040A | Blue  | 5 Slide Trays | 8.9 x 13.4 x 4.9 | 22.6 x 34 x 12.5 | 1/ea |

*\*5 Clasps Included (Clear, Blue, Teal, Purple, Pink)*

## Slide Folder Clasps

| Item No. | Color  | L x W x H<br>in | L x W x H<br>cm | UOM             |
|----------|--|-----------------|-----------------|-----------------|
| 120653   | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | 2.2 x 1.7 x 0.6 | 5.6 x 4.4 x 1.6 | 20 ea color/bag |
| 120654   | Clear <input type="checkbox"/>   |                 |                 | 100 ea/bag      |

*\*Assorted colors (Clear, Gray, Blue, Green, Red)*



120653