

Poulsen & Graf GmbH · Karl-Carstens-Str. 10 · 97877 Wertheim, Germany

April 2024

Technical Information

Dear Ladies and Gentlemen,

The FORTUNA OPTIFIX Dispenser BASIC, SAFETY, SAFETY S, SOLVENT is manufactured out of the below listed materials:


- Soda-lime glass tubing
- PTFE
- PTFE with glass fibre 25 %
- PP
- PVDF
- Sapphire balls
- Ceramic balls
- Borosilicate glass valve seats
- Borosilicate glass tubing

There is no CE marking on the FORTUNA OPTIFIX Dispenser and therefore a CE Certificate is not available.

Each dispenser undergoes a quality test in accordance with ISO 8655-5. The quality certificate and the declaration of conformity, which can be found on the back, are enclosed with the dispensers.

Unfortunately, a food certificate is not available.

Best regards



Sabine Harmon
General Manager / Quality Manager