





# true north<sup>®</sup>

cyro slide box

**Store in ultra-low temperatures down to  
-80°C with confidence**

True North<sup>®</sup> Slide Boxes offer a unique solution for storing slides down to -80°C. Durable polycarbonate construction with a rugged stainless steel clasp provides safe storage and transport for 3x1" (25x75mm) slides. The cryo slide box provides a 100 or 25-place capacity with uniform grooved slide slots to separate and prevent contact.

A permanent identification sheet is included in the hinged lid that corresponds to numbering on the temperature-resistant foam lining in the base. These features provide ideal inventory maintenance ability. The boxes are stackable.

Part No.	Capacity	Hinge	L x W x H	Color
<b>HS15988H</b>	100 slides	Nickle	208 x 175 x 34 mm	Purple 
<b>HS120316</b>	25 slides	Friction fit lid	141 x 88 x 35 mm	Blue 



## Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (25 x 75 mm and 1 x 3 inch)






- Molded, rigid boxes securely hold slides in grooved slots
- Durable ABS plastic boxes are available with either cork or foam lining
- Both the lining and the slots are numbered to correspond to the slide inventory sheet



### 25 Place Slide Box

Easy access to slides with friction-fit cover



Part No.	Capacity	Lining	Hinge	L x W x H	Color
<b>HS15989A</b>	25 slides	Cork	Friction fit	141 x 88 x 35 mm	Blue 
<b>HS15989B</b>	25 slides	Cork	Friction fit	141 x 88 x 35 mm	Green 
<b>HS15989C</b>	25 slides	Cork	Friction fit	141 x 88 x 35 mm	Red 
<b>HS15989E</b>	25 slides	Cork	Friction fit	141 x 88 x 35 mm	White <input type="checkbox"/>
<b>HS15989M</b>	25 slides	Foam	Friction fit	141 x 88 x 35 mm	Blue 
<b>HS15989P</b>	25 slides	Foam	Friction fit	141 x 88 x 35 mm	Red 
<b>HS15989R</b>	25 slides	Foam	Friction fit	141 x 88 x 35 mm	White <input type="checkbox"/>

### 50 Place Slide Box

Thumb-latch lock mechanism











Available in cork lining only

Part No.	Capacity	Lining	Hinge	L x W x H	Color
<b>HS15996A</b>	50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	Blue 
<b>HS15996B</b>	50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	Red 
<b>HS15996C</b>	50 slides	Cork	Thumb-latch	209 x 86 x 35 mm	White <input type="checkbox"/>

### 100 Place Slide Box

Closes securely with a rust-resistant, nickel-plated clasp and pin hinges



Part No.	Capacity	Lining	Hinge	L x W x H	Color
<b>HS15994A</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Blue 
<b>HS15994B</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Green 
<b>HS15994C</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Red 
<b>HS15994D</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Yellow 
<b>HS15994E</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	White <input type="checkbox"/>
<b>HS15994F</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Gray 
<b>HS15994G</b>	100 slides	Cork	Nickle clasp	208 x 175 x 34 mm	Black 
<b>HS15994M</b>	100 slides	Foam	Nickle clasp	208 x 175 x 34 mm	Blue 
<b>HS15994P</b>	100 slides	Foam	Nickle clasp	208 x 175 x 34 mm	Red 
<b>HS15994R</b>	100 slides	Foam	Nickle clasp	208 x 175 x 34 mm	White <input type="checkbox"/>

### 100 Place Tall Slide Box

Stackable with 100-place premium slide boxes

Available in cork lining only  
Additional space for larger 3 x 2 in







Part No.	Capacity	Lining	Hinge	L x W x H	Color
<b>HS15979A</b>	100 slides	Cork	Nickle clasp	208 x 175 x 60 mm	Blue 

## Standard Microscope Slide Boxes

Constructed of durable polypropylene

- Slide boxes hold either 25 or 100 standard slides (25 x 75 mm and 1 x 3 inches)
- Hinged lids are held securely closed with a latch lock
- Stacking feet maximize stability
- Slide identification card
- Unlined and stackable



Part No.	Capacity	Hinge	L x W x H	Color
<b>HS15990A</b>	25 slides	Latch lock	141 x 92 x 30 mm	Blue 
<b>HS15990B</b>	25 slides	Latch lock	141 x 92 x 30 mm	Red 
<b>HS15990C</b>	25 slides	Latch lock	141 x 92 x 30 mm	White <input type="checkbox"/>
<b>HS15991A</b>	100 slides	Latch lock	210 x 169 x 37 mm	Blue 
<b>HS15991B</b>	100 slides	Latch lock	210 x 169 x 37 mm	Red 
<b>HS15991C</b>	100 slides	Latch lock	210 x 169 x 37 mm	White <input type="checkbox"/>