

Customizable Tube Racks

5 mL / 15 mL / 17mm Tube Rack AND
10-17mm Tube Rack

Rotate with ease to accommodate your needs

With bench top space at a premium, it's important for research and clinical labs working across a range of applications to efficiently use the space available. Unlike other tube racks, these multi-purpose racks are specifically designed with the flexibility to hold multiple tube sizes firmly in place to assist technicians constantly alternating between sample vessels.

The rack incorporates two moving hinged shelves that can be rotated to an up or down position to accommodate the size of the tube being used. When working with samples, the shelves can be rotated to an angle position to facilitate more comfortable pipetting and allow better visibility. The tube rack wells have a unique teardrop design that enables easier insertion and better retention of tubes over traditional designs.



Two customizable tube rack versions are available:

For Research Labs: The 5mL/15mL/17 mm tube rack comes with TWO tube segments

For Clinical Labs: The 10-17 mm Tube Rack comes with FOUR interchangeable tube segments.

Additional Features:

- Ergonomic handles for ease of transporting the rack
- Molded grid on both sides of the shelves for quick sample identification and retrieval
- **Autoclavable**
- US Patent Pending




Product Features



Shelves can be rotated to an angle position to facilitate more comfortable pipetting and allow better visibility

Ergonomic handles for ease of transporting the rack

Molded grid on both sides of the shelves for quick sample ID and retrieval



Unique tear drop well design allows for better retention of tubes


Industry	Part No.	# of Segments Included	Tube Sizes	Well Dia.	Wells	L x W x H	Pkg.	Color
Research Labs	HS120167	2	5mL, 15mL, 17 mm	17 mm	20	12.16 x 4.21 x 3.03 in (30.9 x 10.7 x 7.7 cm)	1/ea	Purple Natural <input type="checkbox"/>
Research Labs	HS120467	2	5mL, 15mL, 17 mm	17 mm	20	12.16 x 4.21 x 3.03 in (30.9 x 10.7 x 7.7 cm)	1/ea	Green Blue <input type="checkbox"/>
Clinical Labs	HS120498	2	Up to 17 mm	13 mm	20	12.16 x 4.21 x 3.03 in (30.9 x 10.7 x 7.7 cm)	1/ea	Purple Natural <input type="checkbox"/>
		2	Up to 12 mm	17 mm	20	12.16 x 4.21 x 3.03 in (30.9 x 10.7 x 7.7 cm)		Purple Natural <input type="checkbox"/>

Research Labs

Clinical Labs



HS120167



HS120467




HS120498




HS120498

Comes with 2 tube segments



5mL, 15mL, 17 mm

Comes with 4 tube segments



5mL, 15mL, 17 mm

Up to 12 mm