

Blood Collection Tray

Light weight, heavy on organization

Blood Collection Tray

Organize your phlebotomy, I.V. start or catheter insertion supplies with this spacious, robust and easy to clean carry tray.

- Lightweight easy-to-clean abs construction
- Customize the compartments to suit your personal needs with the dividers provided.
- Comfort grip handle tested to hold up to 10 kg
- Comes complete with a 17 mm well 60 place tube rack

Part No.	L x W x H	Handle Height	Tube Rack
HS120261	16.4 x 12.2 x 3.3 in (41.9 x 31.1 x 8.5 cm)	9.4 in (24.1 cm)	17 mm rack



Comes complete with a 17 mm well, 60 place, tube rack



droplet™

Heathrow Scientific®



Sharps Containers

Easy to use and transport

- Provides safe disposal of syringes, needles and other small sharps
- Stackable to reduce storage space before use
- 1.0 L capacity (quart) with a clear lid to view fill level
- Slide access lid for temporary and final closures for safety and security
- Built-in needle removal port, prevents direct contact with sharps
- Compatible with the droplet® blood collection tray, totemax™ and blood collection tray HS120261

Part No.	L x W x H	Pkg.	Color
HS120177	4.5 x 4.5 x 7.5 in (11.5 x 11.5 x 19.2 cm)	18/pk	Red ■
HS120178	4.5 x 4.5 x 7.5 in (11.5 x 11.5 x 19.2 cm)	18/pk	Yellow ■



x-tourn™ tourniquets

Hi-viz orange alert for safer practice

High-quality disposable tourniquets eliminate roping or curling when used.

Latex free, no rubber proteins, resins or fatty acids that can produce hypersensitive reactions.

- Made from and smooth polyisoprene
- Less abrasive than other brands, uniform thickness prevents shredding and breaking

- Disposable

Part No.	L x W	Pkg.	Color
HS120108	18 x 1 in (45.72 x 2.54 cm)	100/pk	Green ■
HS120109	18 x 1 in (45.72 x 2.54 cm)	25/roll	Green ■
HS120110	18 x 1 in (45.72 x 2.54 cm)	250/pk	Orange ■
HS120111	18 x 1 in (45.72 x 2.54 cm)	25/roll	Orange ■



Available in bulk or on a perforated roll