

# Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss

- Select white, black, or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Can be used with microwave moisture analyzers
- **All are FDA/food-grade compliant**
- **Disposable**



Smooth surfaces provide accurate pour-outs



Manufactured and packaged in a clean room



| Standard weighing boats | Regular Material Part No. | Antistatic Material Part No. | Size   | L x W x H         | Qty.   | Color                                     |
|-------------------------|---------------------------|------------------------------|--------|-------------------|--------|-------------------------------------------|
|                         | <b>HS1420A</b>            | <b>HS1420AA</b>              | Small  | 46 x 46 x 8 mm    | 500/pk | White <input type="checkbox"/>            |
|                         | <b>HS1420B</b>            | <b>HS1420BB</b>              | Medium | 85 x 85 x 24 mm   | 500/pk | White <input type="checkbox"/>            |
|                         | <b>HS1420C</b>            | <b>HS1420CC</b>              | Large  | 140 x 140 x 22 mm | 500/pk | White <input type="checkbox"/>            |
|                         | <b>HS1421A</b>            | <b>HS120222</b>              | Small  | 46 x 46 x 8 mm    | 500/pk | Blue <input checked="" type="checkbox"/>  |
|                         | <b>HS1421B</b>            | <b>HS120223</b>              | Medium | 85 x 85 x 24 mm   | 500/pk | Blue <input checked="" type="checkbox"/>  |
|                         | <b>HS1421C</b>            | <b>HS120224</b>              | Large  | 140 x 140 x 22 mm | 500/pk | Blue <input checked="" type="checkbox"/>  |
|                         | <b>HS1423A</b>            | <b>HS1423AA</b>              | Small  | 46 x 46 x 8 mm    | 500/pk | Black <input checked="" type="checkbox"/> |
|                         | <b>HS1423B</b>            | <b>HS1423BB</b>              | Medium | 85 x 85 x 24 mm   | 500/pk | Black <input checked="" type="checkbox"/> |
|                         | <b>HS1423C</b>            | <b>HS1423CC</b>              | Large  | 140 x 140 x 22 mm | 500/pk | Black <input checked="" type="checkbox"/> |

- Boats are shallow with a wide, flat bottom to resist tipping
- Rounded corners simplify transfer



| Hexagonal weighing boats | Regular Material Part No. | Antistatic Material Part No. | Size    | L x W x H        | Qty.   | Color                          |
|--------------------------|---------------------------|------------------------------|---------|------------------|--------|--------------------------------|
|                          | —                         | <b>HS14251A</b>              | Small   | 38 x 25 x 9 mm   | 500/pk | White <input type="checkbox"/> |
|                          | —                         | <b>HS14251B</b>              | Medium  | 64 x 51 x 19 mm  | 500/pk | White <input type="checkbox"/> |
|                          | —                         | <b>HS14251C</b>              | Large   | 102 x 76 x 22 mm | 500/pk | White <input type="checkbox"/> |
|                          | —                         | <b>HS14251D</b>              | X-Large | 127 x 89 x 31 mm | 500/pk | White <input type="checkbox"/> |

- Bend flexible boats at points to form pour spouts
- Flat bottoms and sloping sides provide accurate transfer



| Pour-boat weighing dishes | Regular Material Part No. | Antistatic Material Part No. | Size   | L x W x H         | Qty.   | Color                          |
|---------------------------|---------------------------|------------------------------|--------|-------------------|--------|--------------------------------|
|                           | <b>HS1419A</b>            | <b>HS120225</b>              | Small  | 43 x 58 x 13 mm   | 500/pk | White <input type="checkbox"/> |
|                           | <b>HS1419B</b>            | <b>HS120226</b>              | Medium | 83 x 132 x 26 mm  | 500/pk | White <input type="checkbox"/> |
|                           | <b>HS1419C</b>            | <b>HS120227</b>              | Large  | 108 x 183 x 26 mm | 500/pk | White <input type="checkbox"/> |

- Pour-spout design easily dispenses samples without flexing the dish
- Flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling



| Square weighing boats | Regular Material Part No. | Antistatic Material Part No. | Size   | L x W x H      | Qty.   | Color                          |
|-----------------------|---------------------------|------------------------------|--------|----------------|--------|--------------------------------|
|                       | <b>HS1422</b>             | —                            | Medium | 90 x 90 x 7 mm | 500/pk | White <input type="checkbox"/> |

- Shallow boat with a wide, flat bottom to resist tipping
- Ideal for solid product weighing



| Diamond weighing boats | Regular Material Part No. | Antistatic Material Part No. | Size   | L x W x H        | Qty.   | Color                                     |
|------------------------|---------------------------|------------------------------|--------|------------------|--------|-------------------------------------------|
|                        | <b>HS1424AA</b>           | <b>HS1425AA</b>              | 5 mL   | 55 x 35 x 6 mm   | 500/pk | White <input type="checkbox"/>            |
|                        | <b>HS1424A</b>            | <b>HS1425A</b>               | 25 mL  | 71 x 46 x 14 mm  | 500/pk | White <input type="checkbox"/>            |
|                        | <b>HS1424B</b>            | <b>HS1425B</b>               | 30 mL  | 78 x 56 x 14 mm  | 500/pk | White <input type="checkbox"/>            |
|                        | <b>HS1424C</b>            | <b>HS1425C</b>               | 100 mL | 119 x 90 x 19 mm | 500/pk | White <input type="checkbox"/>            |
|                        | <b>HS1427AA</b>           | <b>HS1426AA</b>              | 5 mL   | 55 x 35 x 6 mm   | 500/pk | Black <input checked="" type="checkbox"/> |
|                        | <b>HS1427A</b>            | <b>HS1426A</b>               | 25 mL  | 71 x 46 x 14 mm  | 500/pk | Black <input checked="" type="checkbox"/> |
|                        | <b>HS1427B</b>            | <b>HS1426B</b>               | 30 mL  | 78 x 56 x 14 mm  | 500/pk | Black <input checked="" type="checkbox"/> |
|                        | <b>HS1427C</b>            | <b>HS1426C</b>               | 100 mL | 119 x 90 x 19 mm | 500/pk | Black <input checked="" type="checkbox"/> |

- Diamond shape that easily flexes for accurate pouring
- Wide, flat bottom to resist tipping

Flexible boats bend to form pour spout