

# Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss

- Select white, black, or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Can be used with microwave moisture analyzers
- **All are FDA/food-grade compliant**
- **Disposable**



Smooth surfaces provide accurate pour-outs



Manufactured and packaged in a clean room



Standard weighing boats	Regular Material Part No.	Antistatic Material Part No.	Size	L x W x H	Qty.	Color
	<b>HS1420A</b>	<b>HS1420AA</b>	Small	46 x 46 x 8 mm	500/pk	White <input type="checkbox"/>
	<b>HS1420B</b>	<b>HS1420BB</b>	Medium	85 x 85 x 24 mm	500/pk	White <input type="checkbox"/>
	<b>HS1420C</b>	<b>HS1420CC</b>	Large	140 x 140 x 22 mm	500/pk	White <input type="checkbox"/>
	<b>HS1421A</b>	<b>HS120222</b>	Small	46 x 46 x 8 mm	500/pk	Blue <input checked="" type="checkbox"/>
	<b>HS1421B</b>	<b>HS120223</b>	Medium	85 x 85 x 24 mm	500/pk	Blue <input checked="" type="checkbox"/>
	<b>HS1421C</b>	<b>HS120224</b>	Large	140 x 140 x 22 mm	500/pk	Blue <input checked="" type="checkbox"/>
	<b>HS1423A</b>	<b>HS1423AA</b>	Small	46 x 46 x 8 mm	500/pk	Black <input checked="" type="checkbox"/>
	<b>HS1423B</b>	<b>HS1423BB</b>	Medium	85 x 85 x 24 mm	500/pk	Black <input checked="" type="checkbox"/>
	<b>HS1423C</b>	<b>HS1423CC</b>	Large	140 x 140 x 22 mm	500/pk	Black <input checked="" type="checkbox"/>

- Boats are shallow with a wide, flat bottom to resist tipping
- Rounded corners simplify transfer



Hexagonal weighing boats	Regular Material Part No.	Antistatic Material Part No.	Size	L x W x H	Qty.	Color
	—	<b>HS14251A</b>	Small	38 x 25 x 9 mm	500/pk	White <input type="checkbox"/>
	—	<b>HS14251B</b>	Medium	64 x 51 x 19 mm	500/pk	White <input type="checkbox"/>
	—	<b>HS14251C</b>	Large	102 x 76 x 22 mm	500/pk	White <input type="checkbox"/>
	—	<b>HS14251D</b>	X-Large	127 x 89 x 31 mm	500/pk	White <input type="checkbox"/>

- Bend flexible boats at points to form pour spouts
- Flat bottoms and sloping sides provide accurate transfer



Pour-boat weighing dishes	Regular Material Part No.	Antistatic Material Part No.	Size	L x W x H	Qty.	Color
	<b>HS1419A</b>	<b>HS120225</b>	Small	43 x 58 x 13 mm	500/pk	White <input type="checkbox"/>
	<b>HS1419B</b>	<b>HS120226</b>	Medium	83 x 132 x 26 mm	500/pk	White <input type="checkbox"/>
	<b>HS1419C</b>	<b>HS120227</b>	Large	108 x 183 x 26 mm	500/pk	White <input type="checkbox"/>

- Pour-spout design easily dispenses samples without flexing the dish
- Flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling



Square weighing boats	Regular Material Part No.	Antistatic Material Part No.	Size	L x W x H	Qty.	Color
	<b>HS1422</b>	—	Medium	90 x 90 x 7 mm	500/pk	White <input type="checkbox"/>

- Shallow boat with a wide, flat bottom to resist tipping
- Ideal for solid product weighing



Diamond weighing boats	Regular Material Part No.	Antistatic Material Part No.	Size	L x W x H	Qty.	Color
	<b>HS1424AA</b>	<b>HS1425AA</b>	5 mL	55 x 35 x 6 mm	500/pk	White <input type="checkbox"/>
	<b>HS1424A</b>	<b>HS1425A</b>	25 mL	71 x 46 x 14 mm	500/pk	White <input type="checkbox"/>
	<b>HS1424B</b>	<b>HS1425B</b>	30 mL	78 x 56 x 14 mm	500/pk	White <input type="checkbox"/>
	<b>HS1424C</b>	<b>HS1425C</b>	100 mL	119 x 90 x 19 mm	500/pk	White <input type="checkbox"/>
	<b>HS1427AA</b>	<b>HS1426AA</b>	5 mL	55 x 35 x 6 mm	500/pk	Black <input checked="" type="checkbox"/>
	<b>HS1427A</b>	<b>HS1426A</b>	25 mL	71 x 46 x 14 mm	500/pk	Black <input checked="" type="checkbox"/>
	<b>HS1427B</b>	<b>HS1426B</b>	30 mL	78 x 56 x 14 mm	500/pk	Black <input checked="" type="checkbox"/>
	<b>HS1427C</b>	<b>HS1426C</b>	100 mL	119 x 90 x 19 mm	500/pk	Black <input checked="" type="checkbox"/>

- Diamond shape that easily flexes for accurate pouring
- Wide, flat bottom to resist tipping

Flexible boats bend to form pour spout