

<b>CORNING</b>	<b>Life Sciences</b>	Document Number: <b>LSR00157</b>	Rev.: <b>2</b>
Description (Class): <b>Reference, Customer Technical Data Sheet</b> <b>60x15 mm Dish</b>			Page: <b>1 of 2</b>

**Applications:**

- Tissue Culture-treated and Corning® Primaria™ Dishes
  - Growth of adherent cell lines
  - Unimpaired observation of dish contents
- Non-treated dishes
  - Culture of bacteria, fungi, or other microorganisms on a solid substrate (agar) Growth of non-adherent cell lines in liquid medium
- In vitro Fertilization Dish (Cat. No. 353852, only)
  - Recommended for use in IVF labs. Lot specific certification of sterility, non-pyrogenicity, and non-embryotoxicity.

**Features:**

- Flat, distortion free optics
- Lids designed for optimal gas exchange
- Sterile
- Non-pyrogenic
- Packaged in peel open medical style bags
- Durable construction for stable dish manipulation
- Tissue Culture and Corning Primaria-treated
  - Reliable vacuum-gas plasma treatment permanently and consistently modifies the growth surface for cell adherence

**STERILITY:**

Product is gamma irradiated

**REGULATORY COMPLIANCE:**

- Falcon® products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.
- 353852 only - Falcon® products are manufactured under the current ISO 9001 and ISO 13485:2003 Standards, the Medical Device Directive (MDD) 93/42/EEC, and the current FDA Quality System Regulation 21 CFR Section 820.

**RESTRICTIONS:**


- Refer to "Thermoplastics Properties Chart," available at:  
[http://catalog2.corning.com/Lifesciences/media/pdf/an\\_DL\\_226\\_Falcon\\_ThermoPlastics\\_Broch.pdf](http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf)

Description (Class):  
Reference, Customer Technical Data Sheet 60x15 mm DishPage:  
2 of 2

## Related Catalog Numbers:

Cat. No.	Packout Description	Surface
353002	20/bag, 500/case	Tissue Culture Treated
353802	20/bag, 200/case	Corning® Primaria™
351007	20/bag, 500/case	Non-treated
353652	20/bag, 500/case	In Vitro Fertilization product, Tissue Culture-Treated

Growth Area	= 21.29 cm <sup>2</sup>
Working Volume	= 6.0-7.0 mL
Total Volume	= 23 mL



LID



DISH

ITEM	CAT. NO.	DESCRIPTION	DISH AND LID KEY DIMENSIONS mm (inches)						Height of 5 Lidded Dishes
			A	B	C	D	E	F	
1	353002	60 x 15 mm Dish, TC Treated 20/bag, 500/case	54.81 (2.16)	63.34 (2.10)	51.44 (2.03)	13.89 (.55)	56.80 (2.31)	9.63 (.38)	75.18 (2.96)
2	353802	60 x 15 mm Dish, Corning Primaria 20/bag, 200/case	54.81 (2.16)	63.34 (2.10)	51.44 (2.03)	13.89 (.55)	56.80 (2.31)	9.63 (.38)	75.18 (2.96)
3	351007	60 x 15 mm Dish, Non-treated 20/bag, 500/case	54.81 (2.16)	63.34 (2.10)	51.44 (2.03)	13.89 (.55)	56.80 (2.31)	9.63 (.38)	75.18 (2.96)
4	353652	60 x 15 mm Dish, TC Treated 20/bag, 500/case	54.81 (2.16)	63.34 (2.10)	51.44 (2.03)	13.89 (.55)	56.80 (2.31)	9.63 (.38)	75.18 (2.96)

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks).  
© 2013 Corning Incorporated