# Falcon® Pipet Controller



The light, ergonomic, motorized pipet controller is designed for use with glass and plastic serological pipets from 0.5 to 100 mL volume ranges. Conveniently positioned switches change operating modes and speeds quickly to handle different liquid volumes and viscosities. Aspirating and dispensing speed is controlled through the finger triggers. When fully charged, the batteries enable up to 8 hours of continuous use. The large LCD display clearly indicates battery status, pipetting mode and speed.

The Falcon Pipet Controller is supplied with a two-position charging stand, a universal battery charger, and three batteries. Two 0.2  $\mu$ m and two 0.45  $\mu$ m hydrophobic polytetrafluoroethylene (PTFE) filters are supplied with the pipet controller. Additional filters are available as standard accessories. The filters, pipet holder and nosepiece are autoclavable.

### **Features**

- Ergonomic design for exceptional comfort
- ▶ LCD display for continuous visualization of speed, modes and battery status
- Easy access to the battery compartment
- Easy maintenance



# © 2013 Corning Incorporated Printed in U.S.A. 10/13 CLS-DL-LH-364 ROW REV1

# Falcon® Pipet Controller Ordering Information

Cat. No.	Description	
357469	Falcon Pipet Controller, standard version with two 0.2 μm filters and two 0.45 μm filters, 2-position charging stand, universal power supply and set of 3 batteries	
Accessories		
357472	Hydrophobic PTFE filter (0.2 μm), set of 5	
357473	Hydrophobic PTFE filter (0.45 μm), set of 5	
357474	Silicone pipet holder	
357486	Battery (set of 3)	



For additional product or technical information, and a complete listing of our International Offices and Distributors, please visit **www.corning.com/lifesciences** or call 800.492.1110. Outside the United States, please call +1.978.442.2200 or contact your local Corning sales office listed below.

## **Corning Incorporated**

Life Sciences
836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

### Worldwide Support Offices ASIA/PACIFIC Australia/New Zealand t 0402-794-347 China t 86 21 2215 2888 f 86 21 6215 2988 India

t 91 124 4604000

f 91 124 4604099

t 81 3-3586 1996			
f 81 3-3586 1291			
Korea			
t 82 2-796-9500			
f 82 2-796-9300			
Singapore			
t 65 6733-6511			
f 65 6861-2913			
Taiwan			
t 886 2-2716-0338			
f 886 2-2516-7500			

Japan

EUROPE	All Other European
France	Countries
: 0800 916 882	t 31 (0) 20 659 60 51
: 0800 918 636	f 31 (0) 20 659 76 73
Germany : 0800 101 1153 : 0800 101 2427 The Netherlands : 31 20 655 79 28 : 31 20 659 76 73 Jnited Kingdom : 0800 376 8660 : 0800 279 1117	LATIN AMERICA Brasil t (55-11) 3089-7419 f (55-11) 3167-0700 Mexico t (52-81) 8158-8400 f (52-81) 8313-8589

CORNING | FALCON CEllgro AXYGEN. PYREX GOSSELIN