



BSE & TSE Statement

December 13th, 2023

Dear Valued Customer,

This statement is issued to verify that TriForest bottles, flasks, vials, and their caps are made solely from synthetic plastic polymers: [Polycarbonate (PC), Polyethylene Terephthalate (PET), Polyethylene Terephthalate Glycol-Modified (PETG), and Polypropylene (PP)].

No materials of animal origin are used in the production of any TriForest products. They are never in contact with any biological materials at any time during production or product handling. These products are free from known detectable levels of Bovine Spongiform Encephalopathy (BSE) and Transmittable Spongiform Encephalopathy (TSE) causing agents.

All of the materials used in production or processing are certified to be Animal Derived Component Free (ADCF).

Sincerely,

Syed Rehan
Vice President
TriForest Enterprises