
Oosafe® Snap-Cap Tube, Non-Treated Surface

Description (volume) (ml)	Prod. No.	Qty/pc	Qty/cs
14ml	OOPW-OT10	10	500

Intended use

General Intended Use: For the examination of specimens, including blood and tissue donations, derived from the human body.

Specific Intended Use: Aid in collection of samples under negative pressure or by using ordinary methods for the purpose of examination of samples after collection in terms of temperature effect on samples.

Product compartments: 1

Tube Volumes:

Tube compartment 1: 14 ml

Safety

Observe the expiry date and the storage condition on the label.


Do only use optically perfect and undamaged plasticware.

Handling

The package of Oosafe® Snap-Cap Tube is breathable material with no VOC release. After opening the sealing, the tube can be directly put into use.

Although it can be filled up to 14 ml, the recommended volume for use is between 0.5 ml to 12 ml. Follow the scale on the tube.

The cap can be loosely put over the tube to keep the gas changing. The cap can also be pushed down thoroughly for tight sealing.



Discard after use.

Storage: Store at room temperature.