



# Biosigma

Dutscher group

## CLEARStable

Reliable and stable microbial preservation system  
**New beads size!**


*"Simplify your workflow with CLEARStable's larger beads and efficient design"*



1D Version



2D Version



- **Long-Term Preservation:** Ensures viability and stability of microorganisms for extended periods
- **Multiple Freezing Options:** Store microorganisms at  $-20^{\circ}\text{C}$ ,  $-80^{\circ}\text{C}$ , or in nitrogen vapour
- **Larger Beads for Better Handling:** Easier manipulation, improving efficiency
- **Simplified Identification:** Caps and beads in various colors for easy sample identification
- **Flexible Options:** Available with or without cryoprotectant, and with 1D or 2D barcoding








CLEARStable's larger beads make handling easier, allowing for quicker freezing and thawing of your samples



CLEARStable offers a reliable, innovative solution for the stable freezing and long-term preservation of microorganisms at multiple temperature settings. Each vial contains large porous beads, with or without cryoprotectant. Choice of 1D or 2D barcodes for quick sample identification and traceability on boxes and vials.

#### ADVANTAGE

- Extended and stable Preservation
- Multiple Freezing Temperatures
- Easy Handling
- Color-Coded identification
- Customizable Options



**Larger Beads for Easier Handling**

***"Optimize the stability and preservation of your microorganisms with CLEARStable!"***

#### CLEARStable 1D barcoded tube 2ml

CAT. NO.	COLOUR	WITH CRYOPROTECTANT	PACKAGING
CL1D80/AZ	Blue	X	80 vials
CL1D80/G	Yellow	X	80 vials
CL1D80/RO	Red	X	80 vials
CL1D80/V	Green	X	80 vials
CL1D80/BI	White	X	80 vials
CL1D80/MIX	Mixed colour	X	80 vials
CLDR1D80/AZ	Blue		80 vials
CLDR1D80/G	Yellow		80 vials
CLDR1D80/RO	Red		80 vials
CLDR1D80/V	Green		80 vials
CLDR1D80/BI	White		80 vials
CLDR1D80/MIX	Mixed colour		80 vials

#### CLEARStable 2D barcoded tube 2ml

CAT. NO.	COLOUR	WITH CRYOPROTECTANT	PACKAGING
CL2D80/AZ	Blue	X	80 vials
CL2D80/G	Yellow	X	80 vials
CL2D80/RO	Red	X	80 vials
CL2D80/V	Green	X	80 vials
CL2D80/BI	White	X	80 vials
CL2D80/MIX	Mixed colour	X	80 vials
CLDR2D80/AZ	Blue		80 vials
CLDR2D80/G	Yellow		80 vials
CLDR2D80/RO	Red		80 vials
CLDR2D80/V	Green		80 vials
CLDR2D80/BI	White		80 vials
CLDR2D80/MIX	Mixed colour		80 vials

#### Rack 1D Only

CAT. NO.	DESCRIPTION	PACKAGING
BSM58061B1D5/AZ	1D cryo BOX 81 places in PC, Blue grid	5 pcs

#### Rack 2D Only

CAT. NO.	DESCRIPTION	PACKAGING
BSM580612D/5AZ	2D Cryobox 81 places in PC, Blue grid	5 pcs



**Biosigma**  
Dutscher group

Biosigma S.p.A Dominique Dutscher Group  
Via Valletta, 6 | 30010, Cona (Venice) Italy  
Phone ++39 0426 302226 | Fax ++39 0426 302228  
export@biosigmaeu.com | [www.biosigma.com](http://www.biosigma.com)

REV. 2024-10