



Biosigma

Dutscher group

CLEARStab: Advanced Microbial Preservation System

New beads size!

"Simplify your workflow with CLEARStab's larger beads and efficient design"



1D Version



2D Version



- **Long-Term Preservation:** Keeps cultures viable for extended periods
- **Larger Beads for Better Handling:** Easier manipulation, improving efficiency
- **Color-Coded Organization:** Caps and beads in various colors for easy identification
- **Flexible Options:** Available with or without cryoprotectant, and with 1D or 2D barcoding





CLEARStab's larger beads make handling easier, allowing for quicker freezing and thawing of your samples.



CLEARStab offers an innovative, reliable solution for freezing and long-term preservation of microbial cultures. Each vial contains large porous beads, available with or without a glycerol-based cryoprotectant, ensuring maximum culture stability. For streamlined tracking and integration, CLEARStab comes with flexible 1D or 2D barcoding options.

ADVANTAGE

- Extended Preservation
- Easy Handling
- Color-Coded Efficiency
- Flexible Options



Larger Beads for Easier Handling

Optimize your microbial culture preservation with CLEARStab. Contact us today!

CLEARStab 1D barcoded tube 2ml

CAT. NO.	COLOUR	WITH CRYOPROTECTANT	PACKAGING
CL1D80/AZ	Blue	X	80 vials
CL1D80/G	Yellow	X	80 vials
CL1D80/RO	Red	X	80 vials
CL1D80/V	Green	X	80 vials
CL1D80/BI	White	X	80 vials
CL1D80/MIX	Mixed colour	X	80 vials
CLDR1D80/AZ	Blue		80 vials
CLDR1D80/G	Yellow		80 vials
CLDR1D80/RO	Red		80 vials
CLDR1D80/V	Green		80 vials
CLDR1D80/BI	White		80 vials
CLDR1D80/MIX	Mixed colour		80 vials

CLEARStab 2D barcoded tube 2ml

CAT. NO.	COLOUR	WITH CRYOPROTECTANT	PACKAGING
CL2D80/AZ	Blue	X	80 vials
CL2D80/G	Yellow	X	80 vials
CL2D80/RO	Red	X	80 vials
CL2D80/V	Green	X	80 vials
CL2D80/BI	White	X	80 vials
CL2D80/MIX	Mixed colour	X	80 vials
CLDR2D80/AZ	Blue		80 vials
CLDR2D80/G	Yellow		80 vials
CLDR2D80/RO	Red		80 vials
CLDR2D80/V	Green		80 vials
CLDR2D80/BI	White		80 vials
CLDR2D80/MIX	Mixed colour		80 vials

Rack 1D Only

CAT. NO.	DESCRIPTION	PACKAGING
BSM58061B1D5/AZ	1D cryo BOX 81 places in PC, Blue grid	5 pcs

Rack 2D Only

CAT. NO.	DESCRIPTION	PACKAGING
BSM580612D/5AZ	2D Cryobox 81 places in PC, Blue grid	5 pcs



Biosigma
Dutscher group

Biosigma S.p.A Dominique Dutscher Group
Via Valletta, 6 | 30010, Cona (Venice) Italy
Phone ++39 0426 302226 | Fax ++39 0426 302228
export@biosigmaeu.com | www.biosigma.com