

# Microlab 610 Basic Single Syringe Dispenser

Specification Sheet (Part# ML610-DIS )



## Product Specifications

Power Rating	24 VDC, 2.5A
Power Requirements	100-240V 1.5A max 50-60 Hz
Volume Left Syringe	Custom Volume Left Syringe
Program Memory	One program stored in memory while power is on
Syringe Drive Mechanism	Stepper motor driven high precision lead screw with encoder

© 2017 Hamilton Company. All Rights Reserved.

All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

# HAMILTON<sup>®</sup>

Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)  
USA: +1-800-648-5950  
Europe: +41 58 610 10 10

#### Hamilton Americas & Pacific Rim

Hamilton Company Inc.  
4970 Energy Way  
Reno, Nevada 89502 U.S.A.  
Tel: +1-800-648-5950  
Fax: +1-775-856-7259  
[sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

#### Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.  
Via Crusch 8  
CH-7402 Bonaduz, Switzerland  
Tel: +41 58 610 10 10  
Fax: +41-58-610-00-10  
[contact@hamilton.ch](mailto:contact@hamilton.ch)

To find a representative in your area, please visit [hamiltoncompany.com/contacts](http://hamiltoncompany.com/contacts)