

**Arvensis 96 Well PCR Plate, Fully Skirted, Low Profile**

- B-Frame PART N° Dutscher 421118, 421119



Agilent AriaMx qPCR System  
Agilent Robocycler Gradient  
Analytic Jena Alpha SC  
Analytic Jena Flexcycler2, T1 Thermocycler, T Gradient, T One, T Advanced, T Professional series  
Analytic Jena T robot 96, SpeedCycler2 SP, SPR  
Apollo ATC-201  
Biometra T Advanced 96, 96 G, 96 S, 96 SG  
Biometra T Gradient 96  
Biometra T One  
Biometra T Professional series  
Biometra Uno  
Bio-Rad (MJ Research) DNA Engine Opticon and Opticon 2  
Bio-Rad (MJ Research) DNA Engine, DNA Engine Tetrad, Tetrad 2, DNA Engine Dyad, Dyad Disciple  
Bio-Rad (MJ Research) PTC 200  
Bio-Rad C1000, C1000 Touch, S1000  
Bio-Rad C1000, C1000 Touch, S1000  
Bio-Rad CFX Connect, CFX96, CFX96 Touch  
Bio-Rad CFX96 Touch Deep Well  
Bio-Rad Chromo4  
Eppendorf Mastercycler ep realplex  
Eppendorf MasterCycler ep, ep gradient, Pro, Pro S, nexus, nexus gradient, nexus SX1, nexus GSX1  
Eppendorf Mastercycler series  
GE Healthcare Amersham MegaBACE 500, MegaBACE 1000 mark 2  
G-Storm GS1  
Idaho Technology LightScanner  
MJ Research (Bio-Rad) DNA Engine Opticon and Opticon 2  
MJ Research (Bio-Rad) DNA Engine, DNA Engine Tetrad, Tetrad 2, DNA Engine Dyad, Dyad Disciple  
MJ Research (Bio-Rad) PTC100, 200, 220, 221, 225, 240  
peqLab peqSTAR 96X  
Techne Quantica  
Techne TC-412  
Techne TC412, TC512, Genius, Genius Quad, Touchgene, Touchgene Gradient, Flexigene  
Techne TC-PLUS  
Techne Touchgene  
VWR Collection UNO96  
VWR peqlab peqSTAR 96X

**NOTE: We have compiled this list based on plate design and cross-referencing compatibility information available from Eppendorf, Bio-Rad, 4titude, and other sources. Each instrument listed may have variations in block design and configuration. Compatibility is not guaranteed; this is meant as a general guide.**

*All instrument brand names and model names listed above are the trademarks of their respective owners.*

Revised February 10, 2022

**Arvensis 96 Well PCR Plate, Semi-Skirted, Low Profile, Roche style -**

*- B-Frame PART N° Dutscher 421120*

Roche LC96, LC480



**NOTE: We have compiled this list based on plate design and cross-referencing compatibility information available from Eppendorf, Bio-Rad, 4titude, and other sources. Each instrument listed may have variations in block design and configuration. Compatibility is not guaranteed; this is meant as a general guide.**

*All Instrument brand names and model names listed above are the trademarks of their respective owners.*

Revised February 10, 2022

DOMINIQUE DUTSCHER SAS

**Arvensis 96 Well PCR Plate, Semi-Skirted, Standard Profile, Universal style**

- B-Frame PART N° Dutscher 421121, 421122



ABI GeneAmp 9600, 9700  
ABI MiniAmp  
ABI ProFlex  
ABI QuantStudio 3, 5, 6, 7, 12 K Real-time PCR System  
ABI SimpliAmp  
ABI Veriti 96-well  
ABI GeneAmp 2700. 2720  
Agilent (Stratagene) Robocycler Gradient  
Agilent Mx3000P, Mx3005P  
Agilent Surecycler 8800  
Analytica Jena Flexcycler2, T1, Tgradient. Tone, Tadvanced, Tprofessional Gradient/XL, Trobot 96, SpeedCycler2 SP, SPR  
Analytica Jena qTOWER3 (G, touch), Toptical  
Apollo ATC-201  
Bioer Gene Explorer  
Bioer Gene Touch 96  
Biometra T Advanced 96, 96 G, 96 S, 96 SG  
Biometra T Gradient 96  
Biometra T One  
Biometra T Professional series  
Biometra Uno  
BIONEER AllInOneCycler  
Bio-Rad C1000, C1000 Touch, S1000, iCycler, MyCycler, T100  
Bio-Rad CFX96 Touch Deepwell  
Bio-Rad iCycler  
Bio-Rad iQ, iQ-5, myiQ PCR System  
Bio-Rad MyCycler  
Bio-Rad MyiQ, iCycler, IQ, IQ 4, IQ 5  
Bio-Rad PTC 200  
Bio-Rad T100  
Corbett Research (Qiagen) Palm Cycler  
Eppendorf Mastercycler ep realplex  
Eppendorf MasterCycler ep, ep gradient, Pro, Pro S, nexus, nexus gradient, nexus SX1, nexus GSX1  
G-Storm GS1  
HiMedia® Prima-96  
MJ Research (Bio-Rad) Chromo4  
MJ Research (Bio-Rad) Personal, PCT100, 200, 220, 221, 225, 240  
peqLab peqSTAR 96X  
Sensoquest Labcycler  
Stratagene Robocycler Gradient  
Takara Dice touch, Gradient  
Techne Prime, PrimeG, Prime Elite, Prime Elite Satellite, PCRmax Alpha cycler 1, 2, 4  
Techne Genius, Genius Quad, Touchgene, Touchgene Gradient, Flexigene  
Techne TC-412, TC-512. TC Plus  
VWR Collection UNO96

**NOTE: We have compiled this list based on plate design and cross-referencing compatibility information available from Eppendorf, Bio-Rad, 4titude, and other sources. Each instrument listed may have variations in block design and configuration. Compatibility is not guaranteed; this is meant as a general guide.**

*All Instrument brand names and model names listed above are the trademarks of their respective owners.*

Revised February 10, 2022

**Arvensis 96 Well PCR Plate, Semi-Skirted, Standard Profile, ABI style**

- B-Frame PART N° Dutscher 421123, 421124



ABI 3100, 3130XL, 3500, 3500XL, 3730, 3730XL  
ABI 7000, 7300, 7500, 7700, 7900 HT  
ABI GeneAmp 2700, 2720, 9600, 9700  
ABI QuantStudio 3, 5, 6, 7, 12K, ViiA7  
ABI Veriti, Proflex, Simpliamp  
Agilent Mx3000P, Mx3005P, Mx4000  
Analytik Jena Flexcycler2, T1 Thermocycler, Tgradient, Tone, Tadvanced, TProfessional Gradient/XL  
Analytik Jena qTOWER3 G, qTOWER3 touch, Toptical  
Analytik Jena Trobot 96, SpeedCycler2 SP, SpeedCycler2 SPR  
Bioer Technologies Gene Touch 96  
Bio-Rad C1000 Touch, S1000, iCycler, MyCycler, T100  
Bio-Rad MyiQ, iCycler IQ, IQ 4, IQ 5  
Corbett Research (Qiagen) Palm Cycler  
Eppendorf MasterCycler ep, ep gradient, ep realplex, Pro, nexus, nexus gradient, nexus SX1, nexus GSX1  
MJ Research (Bio-Rad) PTC100, 200, 220, 221, 225, 240, Chromo4  
peqLab peqSTAR 96X  
Sensoquest Labcycler  
Takara Dice touch, Gradient  
Techne PCRmax Alpha cycler 1, 2, 4  
Techne Prime, PrimeG, Prime Elite, Prime Elite Satellite  
Techne TC412, TC512, Genius, Genius Quad, Touchgene, Touchgene Gradient, Flexigene

**NOTE: We have compiled this list based on plate design and cross-referencing compatibility information available from Eppendorf, Bio-Rad, 4titude, and other sources. Each instrument listed may have variations in block design and configuration. Compatibility is not guaranteed; this is meant as a general guide.**

*All Instrument brand names and model names listed above are the trademarks of their respective owners.*

Revised February 10, 2022

DOMINIQUE DUTSCHER SAS