

# Corning® Centrifuge Tubes

Designed for Optimum Performance

CORNING

Corning offers a full line of disposable, premium quality conical centrifuge tubes manufactured from high quality, ultra-clear resins. These tubes meet the quality and performance standards demanded by researchers.

▶ **High Speed 15 mL and 50 mL tubes**

- Rated up to 17,000 RCF
- 95 kPA (14 psi) pressure-tested (plug seal cap)
- RNase-/DNase-free
- Nonpyrogenic
- Sterile


▶ **Large Volume 250 mL and 500 mL tubes**

- Rated up to 6,000 RCF
- Best choice for cell harvesting and purification
- Nonpyrogenic
- Sterile

▶ **Simplify your workflow with:**

- Color bands visually identify samples
- Hazardous materials pictograms that comply with OSHA and GHS standards





15 mL and 50 mL centrifuge tubes with CentriStar cap



Easy ID color strip tubes



50 mL centrifuge tube with hazardous materials label

## Ordering Information

### Corning® 15 mL Centrifuge Tubes

Cat. No.	Material	Cap Style	Max. RCF	Feature	Qty/Pk	Qty/Cs
430053	PET	Plug seal	3,600	Conical bottom	25/sleeve	500
430055	PET	Plug seal	3,600	Conical bottom	50/rack	500
430052	PP	Plug seal	12,000	Conical bottom	50/rack	500
430766	PP	Plug seal	12,000	Conical bottom	25/sleeve	500
430790	PP	CentriStar™	12,500	Conical bottom	50/rack	500
430791	PP	CentriStar	12,500	Conical bottom	25/sleeve	500

### Corning 50 mL Centrifuge Tubes

430304	PET	Plug seal	3,600	Conical bottom	25/rack	500
430290	PP	Plug seal	15,500	Conical bottom	25/rack	500
430291	PP	Plug seal	15,500	Conical bottom	25/sleeve	500
430828	PP	CentriStar	17,000	Conical bottom	25/rack	500
430829	PP	CentriStar	17,000	Conical bottom	25/sleeve	500
4558	PP	CentriStar	17,000	Conical bottom	25/Universal rack*	300
<b>New</b> 431526	PP	Plug seal	15,500	Easy ID color strip - green, conical bottom	25/Universal rack	300
<b>New</b> 431527	PP	Plug seal	15,500	Easy ID color strip - Orange, conical bottom	25/Universal rack	300
<b>New</b> 431528	PP	Plug seal	15,500	Easy ID color strip - blue, conical bottom	25/Universal rack	300
<b>New</b> 431525	PP	Plug seal	15,500	Hazardous materials label, conical bottom	25/Universal rack	300
430897	PP	Plug Seal	3,000	Self-standing bottom	25/sleeve	500
430921	PP	CentriStar	3,000	Self-standing bottom	25/sleeve	500

\*New innovative universal rack can hold 50 mL and 15 mL tubes securely, allowing researchers to work with and store both size tubes in the same rack, saving bench and storage space.

### 250 mL and 500 mL Premium Quality Centrifuge Tubes and Support Cushions

Cat. No.	Description	Material	Cap Style	Max RCF	Qty/Pk	Qty/Cs
430776	250 mL Tube	PP	Plug	6,000	6	102
430236	250 mL Support cushion*	PEI	N/A	N/A	N/A	6
431123	500 mL Tube	PP	Plug	6,000	6	36
431124	500 mL Support cushion*	PEI	N/A	N/A	n/a	6

\*Support cushions must be used with the 250 mL and 500 mL centrifuge tubes unless the rotor has appropriately shaped V-bottom holders. PP = Polypropylene, PET = Polyethylene Terephthalate, RCF = Relative Centrifugal Force (x g), PEI = Polyetherimide

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, please visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences). You can contact our Scientific Support Team at [ScientificSupportEMEA@corning.com](mailto:ScientificSupportEMEA@corning.com).

### Corning Incorporated Life Sciences Europe

Corning BV  
Fogostraat 12  
1060 LJ Amsterdam  
The Netherlands  
Phone: +31 (0) 20 659 60 51  
Fax: +31 (0) 20 659 76 73  
[cceurnl@corning.com](mailto:cceurnl@corning.com)  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

### Support Offices

#### EUROPE

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

**All Other European Countries**  
t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

**CORNING** | **FALCON**® **AXYGEN** **GOSSELIN**™