



Equi spaced serrations

Fine spark finish



Thread Details:  
Single start, right hand  
2.0mm pitch



Material: Polypropylene  
Colour: Various  
Gen. Tol: 0.25

**UNCONTROLLED COPY**

NOTE 1. ALL WORK IN ACCORD. WITH WHEATON UK LTD STANDARD SPEC. 2. Ref.	REV	DATE	BY	CHK	APP	Details of last Amendment	
<b>CONFIDENTIAL</b> THIS DRAWING IS THE PROPERTY OF WHEATON UK LTD REPRODUCTION OR DISCLOSURE IS PROHIBITED WITHOUT WRITTEN PERMISSION.	DRAWN P.GANDER					Description	
	CHECKED						15.5mm screw cap
	APPROVED <i>[Signature]</i>					A4	
	DATE 24.01.14						REV
ALL DIMENSIONS IN mm		WHEATON®			1		