

# Dimensions & Weight Of Gel Doc Systems

*System Alone*

## 01. Dims

Height: 585 mm  
Width: 380 mm  
Depth: 475 mm

*System Packed*

## 02. Dims

Height: 860 mm  
Width: 910 mm  
Depth: 560 mm

*System Alone*

## 03. Weight

With Pad: 36 Kg  
Without Pad: 30 Kg

*System Packed*

## 04. Weight

With Pad: 40 Kg  
Without Pad: 34 Kg



DOMINIQUE DUTSCHER SAS