

Tubular Dram Vials

WHEATON® Tubular Dram Vials are manufactured from Type I neutral glass and are available in both **clear** and **amber** in a large range of both tall and squat form sizes. Choose amber glass for protection of light sensitive products. Solid top foil lined closures sold separately provide good barrier properties and resistance to hydrocarbons, oils, ketones and alcohol, making the vials a brilliant packaging solution for oils, perfume samples or tablets.



Ordering Information

Tubular Glass Dram Vials (Amber & Clear)

Code	Design	Glass Type	Capacity (ml)	Diameter x Height (mm)	Cap Size (mm)	Pack Quantity
151061	Tall Form	Clear Type 1 Neutral	1.75	11 x 36	10	882
151062	Tall Form	Clear Type 1 Neutral	3.5	12 X 46	10	666
151063	Tall Form	Clear Type 1 Neutral	7	17 x 58	15	399
151064	Tall Form	Clear Type 1 Neutral	10.5	17 x 67	15	399
151065	Tall Form	Clear Type 1 Neutral	14	20 x 72	18	264
151066	Tall Form	Clear Type 1 Neutral	21	20 x 86	18	231
151067	Tall Form	Clear Type 1 Neutral	28	23 x 96	22	190
152061	Squat Form	Clear Type 1 Neutral	7	20 x 43	18	264
152062	Squat Form	Clear Type 1 Neutral	14	23 x 58	22	190
152063	Squat Form	Clear Type 1 Neutral	28	27 x 72	28	140
151068	Tall	Amber Type 1 Neutral	7	17 x 58	15	399
152064	Squat Form	Amber Type 1 Neutral	14	23 x 58	22	190
152065	Squat Form	Amber Type 1 Neutral	28	27 x 72	28	140

Caps

Code	Cap Material	Cap Liner	Cap Size (mm)	Cap Colour	Pack Quantity
110060-882	Polypropylene	Foil	10	Black	882
110060-666	Polypropylene	Foil	10	Black	666
110061-399	Polypropylene	Foil	15	Black	399
110062-264	Polypropylene	Foil	18	Black	264
110063-190	Polypropylene	Foil	22	Black	190
110064-140	Polypropylene	Foil	28	Black	140