



METRIC STD: (mm)

THREAD TURNS	
Min.	1.50

CODE	MIN.	MAX.
D1	12.50	13.00
⊕ D5	9.40	9.80
D6	7.90	8.30
⊕ D7	5.20	5.75
⊕ H1	45.00	47.00
H3	6.50	7.50


FINISH DETAILS
(ALL DIMENSIONS ARE NOMINAL UNLESS STATED)



⊕ GAUGES NEEDED

CAPACITY (ML) (Nominal)	
Fill Volume	3.10 (85%)
Brim Fill	3.70

BASE DETAILS



CONCAVITY	
Min.	Max.
0.20	

UNCONTROLLED COPY

NOTE:
5 Holed Ring Gauge Available

MATERIAL: TYPE I NEUTRAL CLEAR GLASS

<p>NOTE</p> <p>1. ALL WORK IN ACCORD. WITH WHEATON UK LTD STANDARD SPEC.</p> <p>2. Ref.</p>	REV	DATE	BY	CHK	APP	<p>Details of last Amendment</p>
<p>CONFIDENTIAL</p> <p>THIS DRAWING IS THE PROPERTY OF WHEATON UK LTD REPRODUCTION OR DISCLOSURE IS PROHIBITED WITHOUT WRITTEN PERMISSION.</p>	DRAWN P.GANDER		<p>WHEATON®</p>			<p>Description</p> <p>1 DRAM (3.5ML) LABORATORY THREADED VIAL</p>
	CHECKED	<i>Var Vig 24-6-13</i>				
	APPROVED					
	DATE	24.06.13				
ALL DIMENSIONS IN mm			A4	151062	REV 1	