


REV	APP
3	New layout, base rad reduced, thread detail note added.

Thread Detail:  
Pitch 2mm  
2 Full turns of thread  
Minimum



**UNCONTROLLED COPY**

FREEZE DRY BASE

Material: Amber Neutral Glass

NOTE:  
5 Holed Ring Gauge Available

NOTE  
1. ALL WORK IN ACCORD. WITH WHEATON UK LTD STANDARD SPEC.  
2. Ref.

REV	DATE	BY	CHK	APP
3	26.04.10	PG	SK	Various

Details of last Amendment  
New layout, reduced base rad, thread detail note added

**CONFIDENTIAL**  
THIS DRAWING IS THE PROPERTY OF WHEATON UK LTD REPRODUCTION OR DISCLOSURE IS PROHIBITED WITHOUT WRITTEN PERMISSION.

DRAWN	P.GANDER
CHECKED	<i>Reveling 16-11-12</i>
APPROVED	
DATE	16.11.12
ALL DIMENSIONS IN	mm



Description	
5ml Amber Vial	
A4	102030
REV	4