



WHEATON® COMPLETEPAK™

Designed by you. Crafted by us.



DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

WHEATON® COMPLETEPAK™

The WHEATON® COMPLETEPAK™ product line is a range of crimp top vials, rubber stoppers and seals that are ready-to-use and are off-the-shelf, which can be customized for any application. This reduces our customer's supply chain due to our one source solution capability. Each kit comes with a specific United States Pharmacopeia (USP) & European Pharmacopoeia (EP) certificate showing that the product meets or exceeds critical USP / EP standards that are enforced by the Federal Drug Administration (FDA) & European Directorate for the Quality of Medicines and Healthcare (EDQM). This kind of solution allows for our customers to focus on their core competency which is manufacturing, compounding and developing drug products.

TARGET APPLICATIONS

- Compounding pharmaceuticals
- Aseptic filling of injectable drug products
- Drug discovery
- Drug research
- Freeze drying / lyophilized drug products
- Large and small molecule drug packaging

KEY BENEFITS

- Sterilized, off-the-shelf, ready-to-use (RtU) products
- Variety of Kitted and stand-alone sterile components — a solution for any application
- Kitted solutions feature equal amounts of RtU components for reduced scrap rates
- Complete manufacturing & certificate of analysis traceability
- USP <788>/EP 2.9.19 Particulate Matter in Injections, USP <85>/EP 2.6.1, The Bacterial Endotoxins Test and USP <71>/EP 2.6.14 The Sterility Test, USP <660>/EP 3.2.1 Glass Containers for Pharmaceutical Use, and USP <381>/EP 3.2.9 Elastomeric Closures used in Parenteral Preparations
- COMPLETEPAK™ Amber for light sensitive content packaging
- COMPLETEPAK™ ModPAK Bulk configurations
- Ready-to-Use, Just-in-Time, minimum 180 day shelf life from time of shipment
- Container Closure Integrity (CCI) certification available upon request

MATERIALS

- Vials are designed for parenteral applications and manufactured from low extractable borosilicate glass
- All stoppers are manufactured from an ultrapure bromobutyl latex free formulation with extremely high chemical purity and low gas permeability
- Omniflex stoppers offer low levels of residual binding to maximize extraction of active pharmaceutical ingredients
- Omniflex stopper elasticity delivers a more effective barrier to optimize container closure integrity and content protection
- Igloo lyophilized stoppers are suited for biologics and content with sensitivities to silicone

ENHANCING THE EFFICIENCY OF DRUG DISCOVERY AND DRUG COMPOUNDING



TOP REASONS WHY YOU SHOULD CHOOSE WHEATON® COMPLETEPAK™ FOR YOUR APPLICATION:

1. Products supplied ready-to-use to improve workflow
2. Single source options reduce lead times & simplify administrative effort
3. Focuses your efforts on critical core competencies
4. Eliminates capital investments in sterile processing equipment
5. Frees up valuable clean room space
6. Validated processes & Certificate of Analysis traceability assures quality production outcomes
7. Saves valuable time and money while maintaining a high quality product and streamlining process

WHEATON® COMPLETEPAK™ Kits

Clear Kit Cat. No.	Amber Kit Cat. No.	Qty/ Case	Description (pieces/kit) <i>All pieces packed sterilized</i>	ModPAK Cat. No. (Clear)	ModPAK Cat. No. (Amber)	Cases / ModPAK
W991000-CP	W991000-CPA	1	2mL Vial (278), 13mm Ultra Pure Serum Stopper (285), Blue Seal (285)	W991000-CPMP	W991000A-CPMP	6
W991001-CP	W991001-CPA	1	2mL Vial (278), 13mm Omni Flex Serum Stopper (285), Blue Seal (285)	W991001-CPMP	W991001A-CPMP	6
W991002-CP	W991002-CPA	1	2mL Vial (278), 13mm OmniFlex Igloo Lyo Stopper (285), Blue Seal (285)	W991002-CPMP	W991002A-CPMP	6
W991003-CP	W991003-CPA	1	2mL Vial (278), 13mm Ultra Pure Serum Stopper (285), Red Seal (285)	W991003-CPMP	W991003A-CPMP	6
W991004-CP	W991004-CPA	1	2mL Vial (278), 13mm Omni Flex Serum Stopper (285), Red Seal (285)	W991004-CPMP	W991004A-CPMP	6
W991005-CP	W991005-CPA	1	2mL Vial (278), 13mm OmniFlex Igloo Lyo Stopper (285), Red Seal (285)	W991005-CPMP	W991005A-CPMP	6
W991006-CP	W991006-CPA	1	5mL Vial (120), 20mm Ultra Pure Serum Stopper (125), Blue Seal (125)	W991006-CPMP	W991006A-CPMP	6
W991007-CP	W991007-CPA	1	5mL Vial (120), 20mm OmniFlex Serum Stopper (125), Blue Seal (125)	W991007-CPMP	W991007A-CPMP	6
W991008-CP	W991008-CPA	1	5mL Vial (120), 20mm OmniFlex Igloo Lyo Stopper (125), Blue Seal (125)	W991008-CPMP	W991008A-CPMP	6
W991009-CP	W991009-CPA	1	5mL Vial (120), 20mm Ultra Pure Serum Stopper (125), Red Seal (125)	W991009-CPMP	W991009A-CPMP	6
W991010-CP	W991010-CPA	1	5mL Vial (120), 20mm OmniFlex Serum Stopper (125), Red Seal (125)	W991010-CPMP	W991010A-CPMP	6
W991011-CP	W991011-CPA	1	5mL Vial (120), 20mm OmniFlex Igloo Lyo Stopper (125), Red Seal (125)	W991011-CPMP	W991011A-CPMP	6
W991012-CP	W991012-CPA	1	10mL Vial (104), 20mm Ultra Pure Serum Stopper (125), Blue Seal (125)	W991012-CPMP	W991012A-CPMP	6
W991013-CP	W991013-CPA	1	10mL Vial (104), 20mm OmniFlex Serum Stopper (125), Blue Seal (125)	W991013-CPMP	W991013A-CPMP	6
W991014-CP	W991014-CPA	1	10mL Vial (104), 20mm OmniFlex Igloo Lyo Stopper (125), Blue Seal (125)	W991014-CPMP	W991014A-CPMP	6
W991015-CP	W991015-CPA	1	10mL Vial (104), 20mm Ultra Pure Serum Stopper (125), Red Seal (125)	W991015-CPMP	W991015A-CPMP	6
W991016-CP	W991016-CPA	1	10mL Vial (104), 20mm OmniFlex Serum Stopper (125), Red Seal (125)	W991016-CPMP	W991016A-CPMP	6
W991017-CP	W991017-CPA	1	10mL Vial (104), 20mm OmniFlex Igloo Lyo Stopper (125), Red Seal (125)	W991017-CPMP	W991017A-CPMP	6

WHEATON® COMPLETEPAK™ Sterile Components*

Cat. No.	Description	Qty / Case
Clear Vials		
W987000-S	2mL Clear Sterile Serum Vial	278
W987001-S	5mL Clear Sterile Serum Vial	120
W987002-S	10mL Clear Sterile Serum Vial	104
Amber Vials		
W987000A-S	2mL Amber Sterile Serum Vial	278
W987001A-S	5mL Amber Sterile Serum Vial	120
W987002A-S	10mL Amber Sterile Serum Vial	104

Cat. No.	Description	Qty / Case
Stoppers		
W988003-S	13mm OmniFlex Sterile Serum Stopper	285
W988004-S	13mm Ultra Pure Sterile Serum Stopper	285
W988005-S	13mm OmniFlex Sterile Igloo Lyo Stopper	285
W988006-S	20mm OmniFlex Sterile Serum Stopper	125
W988007-S	20mm Ultra Pure Sterile Serum Stopper	125
W988008-S	20mm OmniFlex Sterile Igloo Lyo Stopper	125
Seals		
W988009-S	13mm Sterile Red Flip Seal	285
W988010-S	13mm Sterile Blue Flip Seal	285
W988011-S	20mm Sterile Red Flip Seal	125
W988012-S	20mm Sterile Blue Flip Seal	125

*Sterile components are available in quantities of 2 or less. For quantities of 3 or more, contact us directly.



Excellence in your hands

DOMINIQUE PUTSCHER SAS

Contact us at [DWK.com](https://www.dwk.com)

© 2021 DWK Life Sciences • The trademarks used are owned by DWK Life Sciences GmbH, DWK Life Sciences LLC or DWK Life Sciences Ltd.