ichemsolute® sorb



Universal Absorbent for oil and chemicals

Th.Geyer

themsolute sorb - Universal Absorbent for oil and chemicals

Absorbing of spilt chemicals in a safe and reliable way and its professional disposal is important to any laboratory.

The new Absorbent chemsolute sorb is offering through innovative technology intelligent functionality decisive advantages compared to other existing Absorbents:

- Fine granules on polymer/mineral base dust-free
- Applicable with any liquid
- No dangerous reactions with chemicals
- Absorbing up to 75-times of its own weight
- Economical and trouble-free dosing
- Saturation level visible by colour intensity
- Barriere function
- Indicator function colour change: acidic (yellow) / alkaline (red)
- Government-certified according to oil absorber guideline type III R suitable for traffic areas, MPA NRW 22 000 8388 10

them solute Sorb is available in the following sizes:





Dosing container 1800 ml

The solution for laboratory benches, fume cupboards or safety cabinets, with two openings for economical sprinkling and forced pouring



Wide-mouth container 4000 ml

More volume - large opening, ideal for use against larger volumes of spilt liquids, simple and economical to store in laboratory cupboards

Art. No. 2440.0001



Storage container 5 kg

Very compact container, stackable, fitting into laboratory cupboards, ideal for use for large volumes of spilt liquids, e.g. for chemical storage areas

Art. No. 2450.0001

Art. No. 2418.0001



Simple to use

chemsolute sorb already shows with minimal dosage optimum absorption results.



chem^{solute®} **SOT**b should be sprinkled equally and extensive on surface of spilt chemicals or oil.



chemsolute sorb provides safe and quick absorption of hazardous substances for clean disposal.



themsolute SOrb embeds the absorbed chemicals into its gel structure and resulting in an almost dry surface.



them^{solute®} **sorb** is an effective liquid barriere, as dam for securing drains or for the prevention of mixing of chemicals.

Th. Geyer

Germany

Th. Geyer GmbH & Co. KG

Dornierstraße 4 D-71272 Renningen

Phone +49 7159 1637-0 Fax +49 7159 18417

 $\pmb{\text{E-Mail}} \quad \text{renningen@thgeyer.de}$

Th. Geyer - Branch Office West -

Kirchstraße 1 D-53797 Lohmar

Phone +49 2246 10098-0 Fax +49 2246 10098-29 E-Mail lohmar@thgeyer.de

Danmark

Th. Geyer Danmark ApS

Industrivej 6A DK-4000 Roskilde

Phone +45 4630 0030 Fax +45 4630 0069 E-Mail sales@thgeyer.dk

Th. Geyer Hamburg GmbH & Co. KG

Oehleckerring 14 D-22419 Hamburg

Phone +49 40 533214-0 Fax +49 40 5321004 E-Mail hamburg@thgeyer.de

Th. Geyer Berlin GmbH

Dillenburger Straße 53 D-14199 Berlin

Phone +49 30 391100-1 Fax +49 303932069 E-Mail berlin@thgeyer.de