

BATCH CERTIFICATE // CERTIFICATE OF ANALYSIS

Description	Microscope slides, ground edges, corners 45°, frosted area white, 76x26x1 mm, pack of 50, LABSOLUTE®
Item no.	7.695 012
Batch no.	all batches
Test report	all test reports

We confirm that all LABSOLUTE® consumables for microscopy (slides & coverslips) comply with the manufacturer's maximum permissible values. Therefore, no errors in use are to be expected that can be attributed to the consumables.

Every batch is tested according to precisely defined parameters by means of visual inspection, testing of physical dimensions, appearance, packaging and function.

The following parameters were tested and met the manufacturer's strict quality requirements:

- Dimensions (length, width and thickness)
- Appearance (colorless, flat, no adhesion, no cloudiness, no pits, grains, stones, scratches or ruptures)
- Package (quantity per box, shapes, labelling and marking of the box)
- Function (refractive index ($n=1,53$), wettability, chemical resistance against air corrosion, boiling test (no cloudiness, bubbles or crackles after boiling), test of the frosting area)

Renningen, 2024-12-18



On behalf of Th. Geyer
Marc Bietz
Product Manager Private Label Product