


NOEX®

TUBES AND SWABS




PLUGS

PRODUCT
PLUGS
Cat. no. ASEPTIC, PL1065PE

DESCRIPTION
Dimensions: fi 16/13,5mm, h=15mm
Material: POLIETYLEN
Weight: 0,6g

PACKAGING
Quantity per PE bag 250 pcs + CARTON FSC
Quantity per carton: 12.000pcs
Weight of carton: 7,5kg




UNIVERSAL PS/TPE SWAB

PRODUCT
UNIVERSAL PS/TPE SWAB
Cat. no. STERIL-SWABU

DESCRIPTION
Dimensions: fi 2,5mm, h=150mm
Material: POLYSTYREN/TPE
Weight: 0,92g

PACKAGING
Quantity per PE bag 250 pcs + CARTON FSC
Quantity per carton: 2.000pcs
Weight of carton: 3,7kg




PHARYNGEL PS/TPE SWAB

PRODUCT
PHARYNGEL PS/TPE SWAB
Cat. no. STERIL-SWABP

DESCRIPTION
Dimensions: fi 2,5mm, h=150mm
Material: POLYSTYREN/TPE
Weight: 0,72g

PACKAGING
Quantity per PE bag 250 pcs + CARTON FSC
Quantity per carton: 2.000pcs
Weight of carton: 3,0kg



NASOPHARYNGEAL PS/TPE SWAB

PRODUCT
NASOPHARYNGEAL PS/TPE SWAB
Cat. no. STERIL-SWABN

DESCRIPTION
Dimensions: fi 4,5mm, h=150mm
Material: POLYSTYREN/TPE
Weight: 0,99g

PACKAGING
Quantity per PE bag 250 pcs + CARTON FSC
Quantity per carton: 2.000pcs
Weight of carton: 3,0kg




TOOLS WORKSHOP



INJECTION MOULDING



DESIGN



SINCE 1986


DESIGN, INJECTION & TOOLING
NOEX®

DESIGN, INJECTION & TOOLING


LABORATORY PRODUCTS




CONICAL CONTAINERS



PETRI DISHES



INOCULATING LOOPS AND SPREADERS



TUBES AND SWABS



ISO 13485




NOEX Sp. z o.o. Sp. J.

ul. Fabianowska 128
62-052 Komorniki, POLAND

+48 61 810 88 15/21

sales@noex.com.pl
www.noex.com.pl

Who we are • Business profile • Our machinery park • How we can help ...



Injection moulding



Tools workshop



Advertising products



Laboratory products



Cooperation production



Technology

NOEX® CONICAL CONTAINERS AND SPOONS




CONICAL CONTAINER 125 ml

PRODUCT
CONICAL CONTAINER 125 ml
Cat. no. ASEPTIC unassembled: PC121/CRG/CVE...
Cat. no. ASEPTIC assembled: PC122PFS
Cat. no. STERILE assembled collectively packing: PC122PP
Cat. no. STERILE assembled individually packing: PC124PI

DESCRIPTION
Dimensions: \varnothing 60mm, h=75mm
Material: POLIPROPYLEN
Weight: 14g
Capacity: 125ml

PACKAGING
Collective: PE BAG + CARTON FSC
Quantity per carton: 200/500szt
Weight of carton: 3,47/8,9g



CONICAL CONTAINER 60 ml

PRODUCT
CONICAL CONTAINER 60 ml
Cat. no. ASEPTIC unassembled: PC60/CVE...
Cat. no. ASEPTIC assembled: PC60PFS
Cat. no. STERILE assembled collectively packing: PC60PP
Cat. no. STERILE assembled individually packing: PC60PI

DESCRIPTION
Dimensions: \varnothing 30mm, h=35mm
Material: POLIPROPYLEN
Weight: 7g
Capacity: 60ml

PACKAGING
Collective: PE BAG + CARTON FSC
Quantity per carton: 200/500szt
Weight of carton: 1,73/4,25g



SPOONS

PRODUCT
SPOON
Cat. no. STERILE white: CUIGHVS
Cat. no. STERILE blue: CUIGHVSB

DESCRIPTION
Dimensions: 150mm
Material: POLIPROPYLEN
Weight: 6g
Capacity: 3ml

PACKAGING
Individual: 1pcs
Quantity per carton: 1000pcs
Weight of carton: 7,0kg

NOEX® PETRI DISHES



Ø65/0 RODAC

PRODUCT
PETRI DISH \varnothing 65/0 RODAC
Cat. no. ASEPTIC: BP650S185Q
Cat. no. STERILE: BP650I185Q

DESCRIPTION
Dimensions: \varnothing 65mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 0
Weight: 8,0g

PACKAGING
Quantity per PE bag: 18 pcs + CARTON FSC
Quantity per carton: 918 pcs
Weight of carton: 9,0kg



Ø65/4 RODAC

PRODUCT
PETRI DISH \varnothing 65/4 RODAC
Cat. no. ASEPTIC: BP654S185Q
Cat. no. STERILE: BP654I185Q

DESCRIPTION
Dimensions: \varnothing 65mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 4
Weight: 8,0g

PACKAGING
Quantity per PE bag: 18 pcs + CARTON FSC
Quantity per carton: 918 pcs
Weight of carton: 9,0kg



Ø55/0 H=14.2mm

PRODUCT
PETRI DISH \varnothing 55/0
Cat. no. ASEPTIC: BP550S155Q
Cat. no. STERILE: BP550I155Q

DESCRIPTION
Dimensions: \varnothing 55mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 0
Weight: 5,7g

PACKAGING
Quantity per PE bag: 15 pcs + CARTON FSC
Quantity per carton: 1005szt
Weight of carton: 7,0kg



Ø55/3 H=14.2mm

PRODUCT
PETRI DISH \varnothing 55/3
Cat. no. ASEPTIC: BP553S155Q
Cat. no. STERILE: BP553I155Q

DESCRIPTION
Dimensions: \varnothing 55mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 3
Weight: 5,7g

PACKAGING
Quantity per PE bag: 15 pcs + CARTON FSC
Quantity per carton: 1005szt
Weight of carton: 7,0kg



Ø90/0 H=14.2mm

PRODUCT
PETRI DISH \varnothing 90/0
Cat. no. ASEPTIC: BP900S255Q
Cat. no. STERILE: BP900I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 0
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg




Ø90/3 H=14.2mm

PRODUCT
PETRI DISH \varnothing 90/3
Cat. no. ASEPTIC: BP903S255Q
Cat. no. STERILE: BP903I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 3
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg




Ø90/0 H=16.2mm

PRODUCT
PETRI DISH \varnothing 90/0
Cat. no. ASEPTIC: BH900S255Q
Cat. no. STERILE: BH900I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=16,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 0
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg




Ø90/3 H=16.2mm

PRODUCT
PETRI DISH \varnothing 90/3
Cat. no. ASEPTIC: BH903S255Q
Cat. no. STERILE: BH903I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=16,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 3
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg

NOEX® PETRI DISHES




Ø90/3 H=16.2mm LASER MARK

PRODUCT
PETRI DISH \varnothing 90/3 LASER MARK
Cat. no. ASEPTIC: BP913S255Q
Cat. no. STERILE: BP913I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 3
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg



Ø90/3 H=16.2mm 2 SECTORY

PRODUCT
PETRI DISH \varnothing 90/3 LASER MARK
Cat. no. ASEPTIC: BP923S255Q
Cat. no. STERILE: BP923I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 2
Vents: 3
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg



Ø90/4 H=14.2mm STABILE

PRODUCT
PETRI DISH \varnothing 90/4 STABILE
Cat. no. ASEPTIC: BP904S255Q
Cat. no. STERILE: BP904I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 1
Vents: 4
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg



Ø90/4 H=14.2mm STABILE 2 SECTORY

PRODUCT
PETRI DISH \varnothing 90/4 STABILE
Cat. no. ASEPTIC: BP924S255Q
Cat. no. STERILE: BP924I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 2
Vents: 4
Weight: 14,0g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,3kg




Ø90/4 H=14.2mm STABILE 3 SECTORY

PRODUCT
PETRI DISH \varnothing 90/4 STABILE
Cat. no. ASEPTIC: BP934S255Q
Cat. no. STERILE: BP934I255Q

DESCRIPTION
Dimensions: \varnothing 90mm, h=14,2mm
Material: POLYSTYREN
Compartment: 3
Vents: 4
Weight: 14,2g

PACKAGING
Quantity per PE bag: 25 pcs + CARTON FSC
Quantity per carton: 600 pcs
Weight of carton: 9,4kg




Ø120/0 H=18mm

PRODUCT
PETRI DISH \varnothing 120/0
Cat. no. ASEPTIC: BP1200S165Q
Cat. no. STERILE: BP1200I165Q

DESCRIPTION
Dimensions: \varnothing 120mm, h=18,0mm
Material: POLYSTYREN
Compartment: 1
Vents: 0
Weight: 23g

PACKAGING
Quantity per PE bag: 16 pcs + CARTON FSC
Quantity per carton: 240 pcs
Weight of carton: 6,5kg




Ø120/4 H=18mm

PRODUCT
PETRI DISH \varnothing 120/3
Cat. no. ASEPTIC: BP1203S165Q
Cat. no. STERILE: BP1203I165Q

DESCRIPTION
Dimensions: \varnothing 120mm, h=18,0mm
Material: POLYSTYREN
Compartment: 1
Vents: 3
Weight: 23g

PACKAGING
Quantity per PE bag: 16 pcs + CARTON FSC
Quantity per carton: 240 pcs
Weight of carton: 6,5kg



Ø140/0 H=20mm

PRODUCT
PETRI DISH \varnothing 140/0
Cat. no. ASEPTIC: BP1400S105Q
Cat. no. STERILE: BP1400I105Q

DESCRIPTION
Dimensions: \varnothing 140mm, h=20,0mm
Material: POLYSTYREN
Compartment: 1
Vents: 0
Weight: 41g

PACKAGING
Quantity per PE bag: 10 pcs + CARTON FSC
Quantity per carton: 110 pcs
Weight of carton: 5,1kg




Ø140/3 H=20

PRODUCT
PETRI DISH \varnothing 140/3
Cat. no. ASEPTIC: BP1403S105Q
Cat. no. STERILE: BP1403I105Q

DESCRIPTION
Dimensions: \varnothing 140mm, h=20,0mm
Material: POLYSTYREN
Compartment: 1
Vents: 3
Weight: 41g

PACKAGING
Quantity per PE bag: 10 pcs + CARTON FSC
Quantity per carton: 110 pcs
Weight of carton: 5,1kg




120/4 H=18mm

PRODUCT
PETRI DISH 120 x 120/4
Cat. no. ASEPTIC: BP120120SQ
Cat. no. STERILE: BP120120SQ

DESCRIPTION
Dimensions: 120x120mm, h=18,0mm
Material: POLYSTYREN
Compartment: 1
Vents: 4
Weight: 30g

PACKAGING
Quantity per PE bag: 15 pcs + CARTON FSC
Quantity per carton: 150 pcs
Weight of carton: 5,5kg

NOEX® INOCULATING LOOPS AND SPREADERS




INOCULATING LOOP 1ul

PRODUCT
INOCULATING LOOP 1ul
Cat. no. quantity per bag 1 pcs: OS1-1
Cat. no. quantity per bag 20 pcs: OS1-20

DESCRIPTION
Dimensions: Length 194mm
Nominal capacity: 1ul
Material: POLYSTYREN
Weight: 1,0g

PACKAGING
Individual: 1 pcs (250pcs) or 20 pcs (500pcs)
Collective: PE BAG + CARTON FSC
Quantity per carton: 2000pcs / 7000pcs
Weight of carton: 4,0kg / 8,8kg




INOCULATING LOOP 10ul

PRODUCT
INOCULATING LOOP 10ul
Cat. no. quantity per bag 1 pcs: OS10-1
Cat. no. quantity per bag 20 pcs: OS10-20

DESCRIPTION
Dimensions: Length 194mm
Nominal capacity: 10ul
Material: POLYSTYREN
Weight: 1,1g

PACKAGING
Individual: 1 pcs (250pcs) or 20 pcs (500pcs)
Collective: PE BAG + CARTON FSC
Quantity per carton: 2000pcs / 7000pcs
Weight of carton: 4,0kg / 8,8kg




SPREADER "L"

PRODUCT
SPREADER "L"
Cat. no. quantity per bag 1 pcs: ETAR1
Cat. no. quantity per bag 5 pcs: ETAR5
Cat. no. quantity per bag 25 pcs: ETAR25

DESCRIPTION
Dimensions: 145x158, \varnothing 3,5mm
Material: POLYSTYREN
Weight: 1,8g

PACKAGING
Individual: 1 pcs (250) or 5 pcs (1250) or 25 pcs (6250)
Collective: PE BAG + CARTON FSC
Quantity per carton: 2000pcs / 10000pcs / 50000pcs
Weight of carton: 4,5kg / 5,8kg / 6,0kg




SPREADER "T"

PRODUCT
SPREADER "T"
Cat. no. quantity per bag 1 pcs: ETAR1T
Cat. no. quantity per bag 5 pcs: ETAR5T
Cat. no. quantity per bag 25 pcs: ETAR25T

DESCRIPTION
Dimensions: 145x158, \varnothing 3,5mm
Material: POLYSTYREN
Weight: 1,8g

PACKAGING
Individual: 1 pcs (250) or 5 pcs (1250) or 25 pcs (6250)
Collective: PE BAG + CARTON FSC
Quantity per carton: 2000pcs / 10000pcs / 50000pcs
Weight of carton: 4,5kg / 5,8kg / 6,0kg

NOEX® TUBES AND SWABS



TL 10

PRODUCT
TUBE TL 10
Cat. no. ASEPTIC PS: TL10PSN
Cat. no. STERILE PS: TL10PSI
Cat. no. ASEPTIC PP: TL10PPN

DESCRIPTION
Dimensions: \varnothing 16mm, h=100mm
Material: POLYSTYREN or POLIPROPYLEN
Weight: PS=5,4g, PP=4,8g
Capacity: 10ml

PACKAGING
Quantity per PE bag: 250 pcs + CARTON FSC
Quantity per carton: 1250 pcs
Weight of carton: 6,5kg / 6,0kg




TR 10

PRODUCT
TUBE TR 10
Cat. no. ASEPTIC PS: TR10PSN
Cat. no. ASEPTIC PP: TR10PPN

DESCRIPTION
Dimensions: \varnothing 13mm, h=100mm
Material: POLYSTYREN lub POLIPROPYLEN
Weight: PS=3,7g, PP=2,8g
Capacity: 7ml

PACKAGING
Quantity per PE bag: 400 pcs + CARTON FSC
Quantity per carton: 1.600 pcs
Weight of carton: PS=6,5kg, PP=5,5kg




TL 65

PRODUCT
TUBE TL 65
Cat. no. ASEPTIC PS: TL65PSN
Cat. no. ASEPTIC PP: TL65PPN

DESCRIPTION
Dimensions: \varnothing 16mm, h=65mm
Material: POLYSTYREN lub POLIPROPYLEN
Weight: PS=3,8g, PP=2,6g
Capacity: 3,5ml

PACKAGING
Quantity per PE bag: 500 pcs + CARTON FSC
Quantity per carton: 2.000 pcs
Weight of carton: PS=7,0kg, PP=6,0kg



TL 75

PRODUCT
TUBE TL 75
Cat. no. ASEPTIC PS: TL75PSN

DESCRIPTION
Dimensions: \varnothing 12mm, h=75mm
Material: POLYSTYREN
Weight: 2,6g
Capacity: 3,5ml

PACKAGING
Quantity per PE bag: 500 pcs + CARTON FSC
Quantity per carton: 3.000 pcs
Weight of carton: 8,5kg




TR 80

PRODUCT
TUBE TR 80
Cat. no. ASEPTIC cap separately: TR80PP
Cat. no. ASEPTIC cap assembled: TR80PPN
Cat. no. STERILE cap assembled: TR80PSI
Cat. no. STERILE cap assembled: TR80PSI

DESCRIPTION
Dimensions: \varnothing 13mm, h=80mm
Material: POLIPROPYLEN
Weight: cap assembled 3,8g, cap separately 2,8g
Capacity: 65ml

PACKAGING
Quantity per PE bag: 500 pcs + CARTON FSC
Quantity per carton: 5.000szt, 10.000szt
Weight of carton: 7,8kg, 13,0kg



PLUGS TR80 Ø13

PRODUCT
PLUG PL 80
Cat. no. ASEPTIC: PL80PPN
Cat. no. STERILE: PL80PSI

DESCRIPTION
Dimensions: \varnothing 13mm
Material: POLIPROPYLEN
Weight: 1,0g

PACKAGING
Quantity per PE bag: 500 pcs + CARTON FSC
Quantity per carton: 5.000szt, 10.000szt
Weight of carton: 7,8kg, 13,0kg




TR 100

PRODUCT
TUBE TR 100
Cat. no. ASEPTIC cap separately: TR100PP
Cat. no. ASEPTIC cap assembled: TR100PPN
Cat. no. STERILE cap assembled: TR100PSI
Cat. no. STERILE cap assembled: TR100PSI

DESCRIPTION
Dimensions: \varnothing 16mm, h=100mm
Material: POLIPROPYLEN
Weight: cap assembled 6,2g, cap separately 4,1g
Capacity: 13ml

PACKAGING
Quantity per PE bag: 250 pcs + CARTON FSC
Quantity per carton: 2500pcs, cap separately - 1500pcs
Weight of carton: cap assembled 8,8kg, cap separately 6,3kg



PLUGS TR80 Ø16

PRODUCT
PLUG PL 100
Cat. no. ASEPTIC: PL100PPN
Cat. no. STERILE: PL100PSI

DESCRIPTION
Dimensions: \varnothing 16mm
Material: POLIPROPYLEN
Weight: 2,0g

PACKAGING
Quantity per PE bag: 250 pcs + CARTON FSC
Quantity per carton: 5.000pcs
Weight of carton: 12,8kg