

HIGH DESERT BIOSCIENCE

CERTIFICATE OF COMPLIANCE

Ethidium Bromide Dropper Bottles

CATALOG # REA-1014

LOT #230335

Expiration Date: 1/1/2029

PRODUCT: Ethidium Bromide Dropper Bottles
QUANTITY: 1 each dropper bottle containing 10 mls
APPEARANCE: Clear, cherry colored solution
FORMULATION (1x): 0.625mg/ml solution

TEST:	Specification	
Solubility (1X concentration)	Clear, haze-free	PASS
Concentration	OD₅₅₀ > 1.400	1.54/PASS
Dropper function	Drop forms with minimal compression of bottle	PASS

REPORT REVIEWED BY:

John King, PhD
Date: December 1, 2023

DOCUMENT #	FORM-PC-5410.100-042299-A
-------------------	---------------------------