


MOUNTING MEDIA

For Processing & Embedding



POLY-MOUNT®

A crystal clear medium providing superior optical clarity and preserving most biological dyes with little or no fading when slides are stored in light-tight containers. The viscosity permits bubble-free applications with minimal spreading and excellent adhesion of coverslips. Poly-Mount® is perfect for mounting and long term storage. It is miscible with toluene and xylene.

Refractive Index: 1.484 | Viscosity: 192 cps

AQUA-POLY/MOUNT

Water-soluble, non-fluorescing mounting medium. Formulated for mounting sections from aqueous solutions. Can also be used for frozen sections, fat stains and immunohistochemical stains when aqueous mounting medium is required. Mountant is easily removed by soaking in water or buffer.

Refractive Index: 1.454 - 1.460 | Viscosity: 390 cps

CITRAMOUNT™ MEDIUM

Butyl acetate and acrylic mixture is the first mounting medium developed that allows coverslipping directly from d-limonene based clearants. Superior results are also obtained with xylene. Recommended for manual coverslipping only.

Refractive Index: 1.484 | Viscosity: 90 cps

POLYGLASS COVERSIPPING MEDIUM

Toluene-based liquid cover-slipping medium designed to eliminate the use of plastic or glass coverslips. Non-yellowing and scratch resistant, medium simplifies the protection of prepared slides. Ideal for uneven thick sections that cause air bubbles under a coverslip with routine procedures. Polyglass can be removed easily by soaking in toluene or xylene.

Refractive Index: 1.48 | Viscosity: 149 cps

ADVANTAGES

- Viscosity permits bubble-free applications with minimal spreading
- Safe and efficient aqueous mounting
- Allows coverslipping directly from d-limonene-based clearants
- Ideal for uneven thick sections that normally cause air bubbles under routine procedures

Mounting Media Products

Cat. #	Description
08381	Poly-Mount® (120 ml, 940 ml)
18606	Aqua-Poly/Mount (20 ml, 100 ml, 1 case - 5 x 20ml)
24214	CitraMount™ Medium (100 ml, 500 ml)
22253	PolyGlass Coverslipping Medium (120 ml, 1 case - 6 x 120 ml)

Order today at [Polysciences.com](https://www.polysciences.com)