

Restriction Enzyme Eag I



Cat.# FG-Eagl Size 500 units Conc.

10 units/ul

Store at -20℃

Supplied with: 10X FastGene® Buffer III (FG-REB3)

10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer

Sterile water

Recognition site



For Research Use Only. Not for use in diagnostic procedures.

ISO9001

Source: *Enterobacter agglomerans*

Reaction conditions

1X FastGene® Buffer III, 37℃ 1X FastGene® FastCut Buffer, 37℃

FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

1X FastGene® Buffer III

50 mM Tris-HCl (pH 7.9 at 25°C) 100 mM NaCl 10 mM MgCl $_2$ 100 μ g/ml BSA

Unit definition

One unit is defined as the amount of enzyme required to digest 1 μg of λ DNA in 1 hour at 37°C in a total reaction volume of 50 μ l.

Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

Dilution buffer:

FastGene® Diluent B

Heat Inactivation

Eag I can be inactivated at 65°C for 20 min.

Methylation sensitivity

dam methylation: Not sensitive dcm methylation: Not sensitive CpG methylation: Sensitive

Relative activity in FastGene® Buffers

FastGene® Buffer I: 10%
FastGene® Buffer II: 25%
FastGene® Buffer III: 100%
FastGene® Buffer IV: 10%
FastGene® FastCut Buffer: 100%

Note

It is an isoschizomer of Xma III. Cleavage of mammalian genomic DNA is blocked by CpG methylation. Activity decreases if buffer pH is not between 7.9 and 9.0 at 25°C.

Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 µg	Xμl
10X FastGene® Buffer III	1 X	5 µl
Eag I	10 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 37°C for 1 hr

- Fast protocol

. dot protoco.		
Component	Final Conc.	Volume
Substrate DNA	1 μg	Xμl
10X FastGene® FastCut Buffer	1 X	5 μΙ
Eag I	10 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 37°C for 15 min

 \times We recommend 5-10 units of enzyme per μg DNA and 10-20 units for genomic DNA in a 1 h digest..