

G Fast Gene

Restriction Enzyme Pvu I



Cat.# FG-Pvul Size 500 units Conc. 10 units/µl

Store at -20°C

Supplied with: 10X FastGene® Buffer III (FG-REB3) 10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer Sterile water

Recognition site

For Research Use Only. Not for use in diagnostic procedures.

Dilution buffer:

FastGene® Diluent B

Heat Inactivation

Methylation sensitivity

dam methylation: Not sensitive *dcm* methylation: Not sensitive CpG methylation: sensitive

Prolonged incubation

A minimum amount of enzyme required to digest 1 μg substrate DNA for 16 hr; 0.13 U.

Relative activity in FastGene® Buffers

FastGene®		25%
FastGene®	Buffer II:	75%
FastGene®		100%
FastGene®	Buffer IV:	50%
FastGene®	FastCut Buffer:	100%

Note

Cleavage of mammalian genomic DNA is blocked by CpG methylation.

Source: Proteus vulgaris

Reaction conditions

1X FastGene[®] Buffer III 37℃ 1X FastGene[®] FastCut Buffer, 37℃

FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

1X FastGene® Buffer III

50 mM Tris-HCl (pH 7.9 at 25°C) 100 mM NaCl 10 mM MgCl₂ 100 μg/ml BSA

Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1 μ g bacteriophage λ at 37°C for 1 hr in 50 μ l reaction mixtures.

Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 µg	Xμl
10X FastGene [®] Buffer III	1 X	5 µl
Pvu I	10 unit	1 µl
Sterile water		up to 50 µl

 \rightarrow Incubate at 37°C for 1 hr

Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 µg	Χ μΙ
10X FastGene [®] FastCut Buffer	1 X	5 µl
Pvu I	10 unit	1 µl
Sterile water		up to 50 µl
1 1 1 1 2700 (15 1		

→ Incubate at 37°C for 15 min

% We recommend 5-10 units of enzyme per μg DNA and 10-20 units for genomic DNA in a 1 h digest.