

Restriction Enzyme Sac I



Cat.# FG-Sacl Size 2,000 units Conc. 20 units/µl

Store at -20℃

Supplied with: 10X FastGene® Buffer I (FG-REB1)

10X FastGene® FastCut Buffer (FG-REBHF)

6X DNA Loading Buffer

Sterile water

Recognition site

For Research Use Only. Not for use in diagnostic procedures.

ISO9001

Dilution buffer

FastGene® Diluent A

Heat Inactivation

Sac I can be inactivated at 65°C for 20 min.

Methylation sensitivity

dam methylation: Not sensitive dcm methylation: Not sensitive CpG methylation: Not sensitive

Prolonged incubation

A minimum amount of enzyme required to digest 1 μg substrate DNA for 16 hr; 0.13 U.

Relative activity in FastGene® Buffers

 FastGene® Buffer I:
 100%

 FastGene® Buffer II:
 75%

 FastGene® Buffer III:
 25%

 FastGene® Buffer IV:
 75%

 FastGene® FastCut Buffer:
 100%

Note

High salts over 10 mM inhibit its activity. Use clean DNA for efficient cleavage. It is not affected by *dam, dcm,* or mammalian CpG methylation.

Source: Streptomyces achromogenes

Reaction conditions

1X FastGene® Buffer I, 37°C 1X FastGene® FastCut Buffer, 37°C

FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 min with FastGene® FastCut Buffer.

1X FastGene® Buffer I

10 mM Bis Tris propane-HCl (pH 7.0 at 25°C) 10 mM MgCl₂ 100 μg/ml BSA

Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1 μg bacteriophage λ (Hind III digestion) at 37°C for 1 hr in 50 μ I reaction mixtures.

Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

Standard reaction condition

- Normal protocol

- Normai protocoi		
Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® Buffer I	1 X	5 μΙ
Sac I	20 unit	1 μΙ
Sterile water		up to 50 μl

- → Incubate at 37°C for 1 hr
- Fast protocol

- rast protocor		
Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® FastCut Buffer	1 X	5 μΙ
Sac I	20 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 37°C for 15 min

 \times We recommend 5-10 units of enzyme per μg DNA and 10-20 units for genomic DNA in a 1 h digest.