

Restriction Enzyme Taq I



Cat.# FG-Taql Size 4,000 units Conc. 20 units/ul

Store at -20℃

Supplied with: 10X FastGene® Buffer III (FG-REB3)

10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer Sterile water

Recognition site

For Research Use Only. Not for use in diagnostic procedures.

ISO9001

Source: Thermus aquaticus YT1

Reaction conditions

1X FastGene® Buffer III 65℃ 1X FastGene® FastCut Buffer, 65℃

FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

1X FastGene® Buffer III

50 mM Tris-HCl (pH 7.9 at 25°C) 100 mM NaCl 10 mM MgCl $_2$ 100 μ g/ml BSA

Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1 μg bacteriophage λ at 65°C for 1 hr in 50 μ l reaction mixtures.

Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

Dilution buffer:

FastGene® Diluent B

Heat Inactivation

Tag I can be inactivated at 80°C for 20 min.

Methylation sensitivity

dam methylation: Conditionally sensitive

dcm methylation: Not sensitive CpG methylation: Not sensitive

Prolonged incubation

A minimum amount of enzyme required to digest 1 μg substrate DNA for 16 hr; 0.5 U.

Relative activity in FastGene® Buffers

FastGene® Buffer I: 50%
FastGene® Buffer II: 100%
FastGene® Buffer III: 100%
FastGene® Buffer IV: 100%
FastGene® FastCut Buffer: 100%

Note

Activity is inhibited by dam methylation partially overlapping its recognition sequence. Incubation at 37°C results in 10% activity obtained at 65°C.

Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® Buffer III	1 X	5 μΙ
Taq I	20 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 65°C for 1 hr

- Fast protocol

- rast protocol		
Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® FastCut Buffer	1 X	5 μΙ
Taq I	20 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 65°C for 15 min

 \times We recommend 5-10 units of enzyme per μg DNA and 10-20 units for genomic DNA in a 1 h digest.