

STYLE #	LINER MATERIAL	COATING MATERIAL	GRIP DESIGN	CUFF STYLE	THICKNESS MM	POWD-ER-FREE	AQL	COLOUR	EN SIZE	LENGTH MM	PACKAGING
73-711	N/A	Neoprene	Bisque fingertips	Rolled beaded	0.175	Yes	0.65	Green	6, 6.5, 7, 7.5, 8, 8.5, 9	310	See below
73-721	N/A	Neoprene	Bisque fingertips	Rolled beaded	0.175	Yes	0.65	Green	6, 6.5, 7, 7.5, 8, 8.5, 9	310	See below



The most advanced glove solution in allergy prevention

PRIMARY INDUSTRIES



IDEAL APPLICATIONS

- Aseptic drug manufacturing
- Biotechnologies
- Medical device manufacturing
- Handling of cytostatics

PACKAGING

73-711:

Pair pack: The glove is packaged one pair per poly wallet, sealed in a poly envelope; 10 pairs per sealed polybag; 5 sealed polybags per master bag; 4 master bags per case liner; 200 pairs per carton. All packaging is printed with cleanroom-compatible IPA-resistant ink. The polybags have an easy-tear opening. The gloves are triple-bagged for superior contamination control.

Bulk pack: 20 right and left hand gloves are packed in separate compartments in an easy-to-open polybag with perforated line opening; 1 polybag per vacuum sealed polybag; 2 polybags per vacuum sealed master polybag; 5 master polybags per case liner; 200 pairs per carton. All packaging is printed with cleanroom-compatible IPA-resistant ink. The gloves are triple-bagged for superior contamination control.

DESCRIPTION

- Made of poly-chloroprene, DermaShield® is the most technically advanced sterile glove available today.
- Free of latex proteins and accelerators, it minimises both Type I and Type IV allergies associated with natural rubber.
- While DermaShield® has the comfort and flexibility of natural rubber gloves, at a similarly competitive price, its neoprene composition offers better chemical resistance. And the new beaded cuff prevents the glove from rolling down, providing increased protection. DermaShield®: for optimal product and worker protection.
- The DermaShield® 73-711 is sterile.
- The DermaShield® 73-721 is available for non-sterile controlled environments.

REMARKS

- Complies to EN455 parts 1 and 2



CATEGORY III

Print date: 21-05-2017

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J

Boulevard International 55, 1070 Brussels, Belgium

Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03

http://www.ansell.eu • E-mail info@ansell.eu