

MSDS INFORMATION

The Spore Ampoules are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200. As such, a Safety Information Sheet has been established to provide general knowledge as to health and safety based upon handling, storage and use of the product.

REF SA1-50-04E, SA1-50-05E, SA1-50-06E, OS1-50-06E, AAD-06E

Identification of Product and Company			
Product Name:	Spore Ampoules of <i>Geobacillus stearothermophilus</i> for Monitoring Steam Processes		
Product Configuration:	Self-Contained Biological Indicator		
Intended Use:	Indicator for Verification of Sterilization Efficacy		
NFPA 704 Hazard Rating	Health: 0	Fire: 0	Reactivity: 0
Prepared By:	J. Goodall	Email:	sales@excelsiorscientific.com
Date Created:	March 25 th , 2019	Date Revised:	NA
Company:	Excelsior Scientific Ltd, Excelsior House, 3 Stow Gardens, Wisbech, Cambridgeshire, PE13 2HS, United Kingdom		
	Emergency Telephone No.: +44 330 124 3050		
Product Composition			
Product consists of a hermetically sealed glass ampule of Tryptic Soy Broth (TSB) modified with a pH indicator containing spores of <i>Geobacillus stearothermophilus</i> in known concentrations.			
Product is free of latex and natural dry rubber products.			
First Aid Measures			
Eyes:	Not applicable		
Skin:	No safety concerns associated with skin contact; wash with soap and water as normal hygiene dictates.		
Ingestion:	Not applicable		
Inhalation:	No safety concerns associated with inhalation.		
Fire Fighting Measures			
Suitable Extinguishing	Use water spray, alcohol resistant foam, dry chemical or carbon		

Media:	dioxide
Fire Hazard:	Not flammable
Explosion Hazard:	Product is not explosive
Handling	
Precautions for Safe Handling:	
Handle in accordance with good industrial hygiene and safety procedures. Wash hands and other exposed areas with soap and water before eating, drinking or smoking and again when leaving work.	
Methods for Containment and Cleaning Up	
For Containment:	Absorb and/or contain spill with inert material, then place in suitable container.
Methods for Cleaning Up:	Clean up spills immediately and dispose of waste safely.
Shipping Information	
Hazard Category:	Non-hazardous
Conditions/Restrictions:	None
Harmonized Tariff Code:	3002.90.5090