| PDS No. 6531xx      | DS No. 6531xx PRODUCT DATA SHEET             |              |             |  | Page 1 of 1     |
|---------------------|--|--------------|-------------|--|-----------------|
| Devision 02         | 60 Well Microtest Plate (Terasaki Plate), PS |              |             |  | 6               |
| Revision 03         |  | Greiner Item | -No. 6531xx |  | greiner bio-one |
| Valid for Item-No.: | 653102                                       | 653180       | 653190      |  |                 |

| 1.  | Description / Specification |   |  |
|-----|-----------------------------|---|--|
| 1.1 | Description                 | PS Microtest Plate (Terasaki Plate), 60 Well, conical well profile,       |  |
|     |                             | flat well bottom, alphanumeric well coding, with lid                      |  |
|     |                             | 653102: Microbatch under Oil Applications, standard                       |  |
|     |                             | 653180, -190: HLA / Microtest Plate, physical surface treatment           |  |
| 1.2 | Dimensions                  | See customer drawing  |  |
| 1.3 | Volume per well             | 10 µl   |  |
| 1.4 | Material / Resin            | Plate and lid: PS (Polystyrene), free of heavy metal                      |  |
| 1.5 | Colour                      | Plate and lid: clear  |  |
| 1.6 | Sterilization               | No  |  |
| 1.7 | Quality Control             | - Raw Material-Control: physical and immunological testing                |  |
|     |                             | - Product-Control: testing of attributive and variable characteristics in |  |
|     |                             | accordance with the valid specification                                   |  |
| 1.8 | Other Information           | For single use only   |  |

| 2.  | Features                 |  |
|-----|--------------------------|--|
| 2.1 | Basic features           |  |
| 2.2 | Temperature range        | -20°C to +60°C   |
| 2.3 | Autoclavability          | No   |
| 2.4 | Centrifugation, max. RCF | N/A  |
| 2.5 | Chemical Resistance      | See homepage:<br>https://www.gbo.com/en_INT/know-how-services/download-center.html |
| 2.6 | Shelf life               | 653102: n/a<br>653180, -190: 1 year after month of production                      |
| 2.7 | Other Information        | -  |

| 3.  | Packaging         | 653102               | 653180   | 653190               |
|-----|-------------------|----------------------|--|----------------------|
| 3.1 | Pieces / Bag      | 10 / bag             | 10 / bag   | 120 / folding carton |
| 3.2 | Pieces / Box      | 580                  | 270  | 480                  |
| 3.3 | Lot-No.           | E JJ MM XXX (manufac | E JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |                      |
| 3.4 | Other Information | -                    |  |                      |

| 4. | Other Information |
|----|-------------------|
|    | -                 |

| Data Sheet subject to change without notice! |         |              |           |                 |
|--|---------|--------------|-----------|-----------------|
|  | Data Sh | neet subiect | to change | without notice! |

| Prior Issue | Drawn            | Approved         | Released         | CONFIDENTIAL: Information contained in this  |
|-------------|------------------|------------------|------------------|--|
| Revision    | Date             | Date             | Date             | document or drawing is confidential and<br>proprietory to Greiner Bio-One GmbH. This             |
| 02          | 26 November 2014 | 28 November 2014 | 28 November 2014 | document may not be reproduced for any   |
| Date        | Name             | Name             | Name             | reason without written permission from Greiner<br>Bio-One GmbH. All rights of design, invention, |
| 18.12.2009  | S. Kaelberer     | Dr. U. Honisch   | A. Schulz        | and copyright are reserved.  |

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.