

New Plate Printing

Article No. 657160

CELLSTAR® Multiwell Plate, 6 Well, PS, clear, TC

In addition to the established design and proven quality the above mentioned plate will come in future with a black printing as novel feature.

The printing was designed to improve the identification of the respective plate as well as its usability. It includes:

- **Alphanumeric coding for faster orientation**
- **Type of surface treatment for simple and fast plate identification**
- **Lot number and expiration date**
(For technical reasons, bar coded plates can just be printed with alphanumeric coding.)

Furthermore, a new frosted field at the side of the plate enables individual labelling by the customer.

The existing article number **657160** will be retained for the new plates. For a smooth transition, boxes containing the printed version are labeled with “**New Design**”.



6 well multiwell plate with black printing of:

- Alphanumeric coding
- Lot number
- Expiry date
- Surface treatment

For further information and/or sample ordering please visit our website www.gbo.com or contact us.

Germany (Main office): info@de.gbo.com | Austria: office@at.gbo.com | Belgium: info@be.gbo.com | Brazil: office@br.gbo.com | China: office@cn.gbo.com
France: accueil@gbo.com | Hungary: office@hu.gbo.com | Italy: office@it.gbo.com | Japan: info@jp.gbo.com | Netherlands: info@nl.gbo.com
Portugal: info@vacuette.pt | Spain: info@vacuette.es | UK: info@uk.gbo.com | USA: info@us.gbo.com