




Valid for Item-No.: 664910, 664920, 664940, 664950



Curvature +0.15  
-0.70



$\phi$  89,05 <sup>+0.1</sup>/<sub>-0.2</sub>



Curvature +0.50  
-0.15

$\phi$  87,05

"All dimensions in mm"

Customer drawing subject to change without notice!

Uncontrolled Copy

| Prior Issue        | Drawn              | Approved           | Released           | <small>CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.</small> |
|--------------------|--------------------|--------------------|--------------------|---|
| Revision<br>2      | Date<br>2011/03/07 | Date<br>2015/04/02 | Date<br>2015/04/02 |   |
| Date<br>2012/12/12 | Name<br>forb       | Name<br>forb       | Name<br>sonntag    |   |