

Shorr Solution

Code product: **361100**

Specifications:

- **Appearance:**
- **Absorption Maxima in water:**

Black solution with purplish blue tints

$498 \text{ nm} \leq \lambda_1 \leq 508 \text{ nm}$

$622 \text{ nm} \leq \lambda_2 \leq 626 \text{ nm}$

Applications:

Major Application in Hormonal Cytology: used as a cytoplasmic staining solution in Harris-Shorr staining, progressive and regressive variants, for the bichromic staining of vaginal smears. Also used in Shorr staining for the quick bichromic staining of vaginal smears.

Other Application in Cytology: Used in Haemalum-Shorr staining and in simplified Harris-Shorr staining for spermocytograms (morphological analysis of spermatozoa).

Updated: 25/04/2013