

# Corning® Gosselin™ Round HDPE Storage Bottles

High density polyethylene bottles with wide neck

CORNING

- ▶ HDPE storage bottles are easier to handle and are ideal for sampling and storage.
- ▶ HDPE resin is lightweight and break-resistant
- ▶ Bottles are leak-proof

## Key Features

- ▶ Good for storage of liquid, semi-liquid, and solid forms
- ▶ Available in six sizes: 50 mL, 100 mL, 150 mL, 250 mL, 500 mL, and 1L
- ▶ Large neck opening provides easier access for liquid transfer
- ▶ Sterile
- ▶ Raw materials are food contact compliant



## Ordering Information

### Corning® Gosselin™ Round HDPE Storage Bottles

Cat. No.	Description	Qty/Cs
LR50A-18	50 mL, 27 mm red screw cap, sterile	600
LR100-12	100 mL, 37 mm red screw cap, sterile	335
LR150A-14	150 mL, 37 mm red screw cap with natural seal, sterile	250
LR250-09	250 mL, 37 mm red screw cap, sterile	145
LR500-04	500 mL, 58 mm red screw cap, sterile	140
LR1000-08	1L, 58 mm red screw cap, sterile	68

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# CORNING

**Corning Incorporated**  
Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

#### ASIA/PACIFIC

**Australia/New Zealand**  
t 61 427286832

**China**  
t 86 21 3338 4338  
f 86 21 3338 4300

**India**  
t 91 124 4604000  
f 91 124 4604099

#### Japan

t 81 3-3586 1996  
f 81 3-3586 1291

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6572-9740  
f 65 6735-2913

**Taiwan**  
t 886 2-2716-0338  
f 886 2-2516-7500

#### EUROPE

CSEurope@corning.com

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 020 655 79 28  
f 020 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

#### All Other European Countries

t +31 (0) 206 59 60 51  
f +31 (0) 206 59 76 73

**LATIN AMERICA**  
grupoLA@corning.com

**Brasil**  
t 55 (11) 3089-7400

**Mexico**  
t (52-81) 8158-8400