

The logo for ClearLine is displayed in white text on a red, slanted rectangular background. The word "ClearLine" is written in a sans-serif font, with a registered trademark symbol (®) to the upper right of the "e".

**ClearLine<sup>®</sup>**

## STERILE SALINE SOLUTION

---

### INTENDED USE

Sterile saline solution is an isotonic diluent use for dilutions in food microbiology

### FORMULA

Ingredients in grams per liter of purified water

Sodium chloride	8.50
-----------------	------

### STORAGE

Tubes, bottles and bags: 2 - 25°C

The expiration date on the product label applies to the product in its intact packaging when stored as directed.

### LIMITATION OF THE PROCEDURE

This product is for laboratory use only.

DOMINIQUE DUTSCHER S.P.A.