2mag STIRRING SYSTEMS FOR CELL CULTURES AND TISSUE CULTURES

Designed for careful and protective mixing of cell suspensions and culture broths, developed for highest demands as well as warming-free and uncompromising continuous operations for e.g. in CO₂-incubators.



2mag bioMIXdrive 2/3/4

2mag magnetic (e)motion

Schragenhofstr. 35 80992 München

Tel. +49 (89) 14 33 42 52 Fax +49 (89) 14 33 43 69

magnetic motion

			0	
Туре	bioMIXdrive 2	bioMIXdrive 3	bioMIXdrive 4	bioMIXdrive 6
Stirring points	2	3	4	6
Stirring point distance	140 130			
Stirring volume /point	5 - 5,000 ml			
Speed range	5 - 250 rpm			
Stirring power (max)	6 W / position			
Material housing	Stainless steel			
Operation conditions	Operation conditions -10 up to +50 °C (at 100% humidity)			
Operation conditions	Operation conditions +50 °C submerged in water			
Operation voltage	48 VDC			
Dimensions (WxDxH)	130 x 270 x 38	130 x 410 x 38	270 x 270 x 38	245 x 375 x 38 mm
Protection class	IP 68			
Weight (gross)	approx. 2.5 kg	approx. 3.5 kg	approx. 4.7 kg	ca. 7,8 kg
Order No.	80002	80003	80004	80006

Individual units and dimensions on request

3 years warranty on material and manufacturing. We reserve the right to make technical alterations. We do not assume liability for errors in printing or typeing.