

# Azlon®

New to the  
Azlon® Wash Bottle  
range...

Get in touch  
for your free  
sample!

## The Azlon® Round Integral Wash Bottle



## Introduction

We are proud to introduce the new Round Integral Wash Bottle to the Azlon® range. With this design you can be sure to dispense liquids in a safe and efficient manner with extra protection from contamination.

The new Azlon® Round Integral Wash Bottle displays chemical names and key safety information in four European languages and also meets the latest requirements in European safety legislation – plus all of the regular features and benefits of this popular laboratory essential.

## Key Features

- Wide mouth for easy filling
- Colour-coded lids and contents
- Precision jet tip for fine control
- Choice of indelibly printed solvent labels
- Custom printing available on request

Traditional wash bottles consist of several components – a bottle, a cap and an inner tube. The inner tube can be a source of contamination when refilling as it is typically placed on a lab bench surface.

Chemical names in English, French, German and Spanish eliminate the need for multiple wash bottles across European language speaking countries.

Chemical formulae and CAS number for easy identification of solvent contained.

Globally Harmonised Symbols reflects the latest changes in European legislation including new pictograms.

Signal word (i.e. 'Danger'), to identify the severity of the hazard statements.

The Azlon® Round Integral Wash Bottle has an integral inner straw that eliminates the need to place the inner tube down when refilling, minimising the risk of contamination.

Vented cap with PTFE membrane allows venting and prevents pressure build ups.

Precision jet tip ensures fine control.

Moulded graduations for approximate volume measurement.

Integral inner tube prevents bench-derived contamination when refilling.

## Technical and Ordering Information

### Azlon® Round Integral Wash Bottle

- Wide neck wash bottle with red, yellow, blue, green or orange cap
- 500ml capacity
- Fitted with colour-coded one-piece PP closure that maintains its shape in use
- Designed with integral spout to prevent contamination

| Product Code | Closure Colour                                   | Pack Quantity |
|--------------|--------------------------------------------------|---------------|
| WGW800P      | Mixed pack – red, yellow, blue, green and orange | 5             |

### Azlon® Multilingual Round Integral Wash Bottle

- Multilingual, wide neck wash bottle
- 500ml capacity
- Indelibly printed in English, French, German and Spanish
- Choose from five solvent labels all including the new Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Fitted with colour-coded one-piece PP closure that maintains its shape in use
- Designed with integral spout to prevent contamination

| Product Code | Print Label     | Closure Colour | Pack Quantity |
|--------------|-----------------|----------------|---------------|
| WGW801P      | Distilled Water | White          | 5             |
| WGW802P      | Acetone         | Red            | 5             |
| WGW803P      | Methanol        | Green          | 5             |
| WGW804P      | Isopropanol     | Blue           | 5             |
| WGW805P      | Ethanol         | Orange         | 5             |



For your free sample,  
please contact  
[enquiries@scilabware.com](mailto:enquiries@scilabware.com)  
quoting **AZRIWB18**



Email us with your questions at [enquiries@scilabware.com](mailto:enquiries@scilabware.com)

SciLabware Limited, Unit 4, Riverside 2, Campbell Road, Stoke-on-Trent, ST4 4RJ, UK  
Tel: +44 (0)1782 444406 | Fax: +44 (0)1782 940436



SciLabware is a company of **DWK Life Sciences**