




AC- 70107x	CUSTOMER DRAWING	Page 1 of 1
Revision 01	96 Well ELISA-Plate, PS, 8 x F12 Strips, F-Bottom	 greiner bio-one
	Item-No. 701 07x	



NOT TO SCALE



Cat.-No.	Description	A	B	C	D	E	F	G	H	I
701 070	ELISA Plate	127.8 +0.1/ -0.2	85.6 +0.1/ -0.2	14.4 +0.3/ -0.2	63.12	99.11	109.35 ± 0.1	Ø 7.0 +0.05/ -0.15	10.9 -0.15	12.41 ± 0.1
701 071	ELISA Plate	127,8 +0.1/ -0.2	85.6 +0.1/ -0.2	14.4 +0.3/ -0.2	63.12	99.11	109.35 ± 0.1	Ø 7.0 +0.05/ -0.15	10.9 -0.15	12.41 ± 0.1

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date April 5, 2006	Date April 6, 2006	Date April 6, 2006	
Date	Name F. Falter	Name F. Sonntag	Name F. Sonntag	