

CORNING

Trypsin and Trypsin/EDTA Solutions

Formulation

Cat. No.	25-050	25-051	25-052	25-053	25-054
Description	Liquid, 1x 0.25% Trypsin in HBSS	Liquid, 1x 0.05% Trypsin 0.53 mM EDTA in HBSS	Liquid, 1x 0.05% Trypsin 0.53 mM EDTA in HBSS	Liquid, 1x 0.25% Trypsin 2.21 mM EDTA in HBSS	Liquid, 10x 2.5% Trypsin in HBSS
Units	mg/L	mg/L	mg/L	mg/L	
Components					
D-Glucose	1000.00	1000.00	1000.00	1000.00	1000.00
EDTA • 4Na • 2H ₂ O	--	221.00	221.00	920.00	--
KCl	400.00	400.00	400.00	400.00	400.00
KH ₂ PO ₄	60.00	60.00	60.00	60.00	60.00
NaCl	8000.00	8000.00	8000.00	8000.00	8000.00
NaHCO ₃	350.00	350.00	--	--	350.00
Na ₂ HPO ₄ (anhydrous)	47.70	47.70	47.70	47.70	47.50
Phenol red • Na	10.00	10.00	10.00	10.00	--
Trypsin 1:250	2500.00	500.00	500.00	2500.00	25000.00



At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please visit www.cellgro.com or call 1.800.235.5476.

Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
t 800.235.5476
t 703.471.5955
f 703.467.9851
www.cellgro.com

CORNING | **FALCON** | **cellgro**® | **AXYGEN** | **PYREX**® | **GOSELIN**®

Corning cellgro products are manufactured by Mediatech, Inc.
For a listing of Corning trademarks, visit us at www.corning.com/about-us/trademarks.html
All other trademarks in this document are the property of their respective owners.