

# NEPTUNE®

Tools for Life Sciences

## Optimus Rack Upgrade

### Functional Tabs

- Simply squeeze and pull up for easy one handed snap card removal with reloads.

### Locking Clasp

- Secures tips during transport.
- Clicks to ensure rack is securely closed.




### Double Ball Hinge

- Upgraded hinge design withstands repeatable autoclaves.

### Modified Hole Design

- Prevents rack lift up when removing tips.



### Structural Support

- New SBS base footprint makes the rack more stable when attaching tips, eliminating the need for a Suretray.

### Reusability

- The base is engineered with pyramidal wells and an SBS footprint making it reusable as a 96-well automation reservoir.

## See what researchers and scientists have said about the Optimus Rack upgrades

*"It is great! Very user-friendly."*

*"The hinged design works well with multichannel pipets and the design of the box keeps the tips from all falling out when you place them on your pipette."*

*"I love that it is eco-friendly! Helps lab be greener, BIG positive!!"*

*"The locking clasp is very sturdy, superior to other brands."*

[www.neptunescientific.com/optimus](http://www.neptunescientific.com/optimus)