

Frickenhausen, 4 November 2020

INFOMAILING


PRODUCT MODIFICATION: FRAMES OF 96 WELL ELISA STRIP PLATES

Dear Sir or Madam,


We are continuously striving to optimise our product portfolio and to further improve existing products. Therefore we would like to inform you about a minor design change for the frame of our 96 well ELISA strip plates. The modified frames are expected to be available **from late Q4/2020**.

We will introduce **4 slits to the inner corners of the frame**. This helps **minimise the abrasion during the assembly** of frame and strips.

CURRENT FRAME DESIGN



NEW FRAME DESIGN



Items affected by this change:

Item No.	Product Description
754061	STRIP PLATE, 6xU16, PS, U-BOTTOM, CLEAR, MICROLON®, HIGH BINDING, 10 PCS./BAG
754070	STRIP PLATE, 6xU16, PS, U-BOTTOM, CLEAR, MICROLON®, MED. BINDING, 10 PCS./BAG
754074	STRIP PLATE, 6xU16, U-BOTTOM, WHITE, HIGH BINDING
754101	PS-FRAME, WHITE 128.0/86 MM
756070	STRIP PLATE, 6xF16, PS, F-BOTTOM, CLEAR, MICROLON®, MED. BINDING, 10 PCS./BAG
756071(-XXX)	STRIP PLATE, 6xF16, PS, F-BOTTOM, CLEAR, MICROLON®, HIGH BINDING
762070(-XXX)	STRIP PLATE, 12xF8, PS, F-BOTTOM, CLEAR, MICROLON®, MED. BINDING
762071(-XXX)	STRIP PLATE, 12xF8, PS, F-BOTTOM, CLEAR, MICROLON®, HIGH BINDING
762074	STRIP PLATE, 12xF8, PS, F-BOTTOM, WHITE, LUMITRAC, HIGH BINDING, 10 PCS./BAG
762075	STRIP PLATE, 12xF8, PS, F-BOTTOM, WHITE, LUMITRAC, MED. BINDING, 10 PCS./BAG
762076	STRIP PLATE, 12xF8, PS, F-BOTTOM, BLACK, FLUOTRAC, MED. BINDING, 10 PCS./BAG
762077	STRIP PLATE, 12xF8, PS, F-BOTTOM, BLACK, FLUOTRAC, HIGH BINDING, 10 PCS./BAG
767070	STRIP PLATE, 12xU8, PS, U-BOTTOM, CLEAR, MICROLON®, MED. BINDING, 10 PCS./BAG
767071	STRIP PLATE, 12xU8, PS, U-BOTTOM, CLEAR, MICROLON®, HIGH BINDING, 10 PCS./BAG

If you have questions, please do not hesitate to contact your export manager, local subsidiary or GBO distributor.

Kind regards,



i.V. Michael Wirth
 Director Marketing & Product Management